



# RDM Industrial Products Inc.

QUALITY • SERVICE • TIMING - phone (877) 777-9130



Bench Top Laminar Flow Unit



Standard 72" Laminar Flow Workstation



Custom...

Laminar Flow Hood



Vertical Hood





# RDM Industrial Products Inc.

QUALITY • SERVICE • TIMING - phone (877) 777-9130



## Model LH-102P Horizontal Clean Air Workstation Hood:


**This Horizontal Laminar Flow Unit is Great for Operations that Do Not require overhead clearance...**

- Standard Width Sizes include 36", 48", 60", 72" and 96"
- Customizable (RDM can accommodate custom needs)
- Light & Blower Switch are combined into one so that user knows when unit is in operating mode.
- Fluorescent Light Fixture(s)
- Plugs into standard 115v Outlet
- Elec. Junction Box with 12' cord
- Optional Features (Not Included):
  - Full-Height, Vertical Slat Curtain available
  - ULPA Filter(s) available
  - ESD Work Surface
  - Static Control Ionizing Grid(s) available
  - Power Outlet Strip(s) – 8 & 16 Outlet



## Sample Products

### LFH-102P-6 HORIZONTAL LAMINAR FLOW HEPA WORKSTATION



TOTAL UNIT QTY.: \_\_\_\_\_  
OVERALL DIMENSIONS: 3'-9"D x 6'-5"W x 5'-1 1/2"H

#### Dimensions:

Dimension A (WIDTH) -	6'-5"
Dimension A1 (TOP WIDTH) -	6'-0 7/8"
Dimension B (DEPTH) -	3'-9"
Dimension B1 (TOP DEPTH) -	2'-3 7/8"
Dimension C (HEIGHT) -	5'-1 1/2"
Dimension D (TOP HEIGHT) -	3'-0"
Dimension E (CLR. OPEN'G) -	1'-10 1/2"

#### Standard Features:

- Dual HEPA Fan Filter Units w/ Speed Control
- Durable Welded Tube Steel Construction
- Plastic-Laminated Main Work Surface (Shelf Optional)
- White Polyurethane Paint
- 1/2" Side Panels
- Appx. 400# with Laminate Top
- 12' Long, 115v Power Cord

#### Surface Options:

- Stainless Steel, Type 304, #4 Finish
- Perforated Stainless Steel, Type 304, #4 Finish
- ESD Laminate Color \_\_\_\_\_
- 180 Degree Rolled Front Edge
- Power Strip  6  8  15
- Clean Room Compatible


#### General Options:

- Cleanroom Shroud
- Casters/Rollers  W/ Brakes
- Cleanroom Casters  W/ Brakes
- ABS Plastic Modular Drawer System
- Metal File Drawers  Metal Utility Drawers
- Full Front Panel
- Lower Shelf  HALF  FULL

**RDM** RDM INDUSTRIAL PRODUCTS, INC.  
1652 Watson Court  
Milpitas, California 95035  
Ph:(408) 945-8400 Fax(408) 945-8433

## Sample Products

### LFH-102P-6 HORIZONTAL LAMINAR FLOW HEPA WORKSTATION



TOTAL UNIT QTY.: \_\_\_\_\_  
**OVERALL DIMENSIONS: 3'-9"D x 6'-5"W x 5'-10 1/2"H**

Light & HEPA FFU are both listed.  
 Total Power Draw appx. 5amps  
 (operating) 9amps (start-up)

#### Dimensions:

Dimension A (WIDTH) -	6'-5"
Dimension A1 (TOP WIDTH) -	6'-0 7/8"
Dimension B (DEPTH) -	3'-9"
Dimension B1 (TOP DEPTH) -	2'-3 7/8"
Dimension C (HEIGHT) -	5'-10 1/2"
Dimension D (TOP HEIGHT) -	3'-0"
Dimension E (CLR. OPEN'G) -	1'-10 1/2"

#### Standard Features:

- Dual HEPA Fan Filter Units w/ Speed Control
- Durable Welded Tube Steel Construction
- Plastic-Laminated Main Work Surface (Shelf Optional)
- White Polyurethane Paint
- 1/2" Side Panels
- Appx. 400# with Laminate Top
- 12' Long, 115v Power Cord

#### Surface Options:

- Stainless Steel, Type 304, #4 Finish
- Perforated Stainless Steel, Type 304, #4 Finish
- ESD Laminate Color \_\_\_\_\_
- 180 Degree Rolled Front Edge
- Power Strip  6  8  15
- Clean Room Compatible

#### General Options:


- Cleanroom Shroud
- Casters/Rollers  W/ Brakes
- Cleanroom Casters  W/ Brakes
- ABS Plastic Modular Drawer System
- Metal File Drawers  Metal Utility Drawers
- Full Front Panel
- Lower Shelf  HALF  FULL

RDM

**RDM INDUSTRIAL PRODUCTS, INC.**  
 1652 Watson Court  
 Milpitas, California 95035  
 Ph:(408) 945-8400 Fax(408) 945-8433

## Sample Products

### LFH-102P-8 HORIZONTAL LAMINAR FLOW HEPA WORKSTATION



TOTAL UNIT QTY.: \_\_\_\_  
OVERALL DIMENSIONS: 3'-9"D x 8'-5"W x 5'-1 1/2"H

Light & HEPA FFU are both listed.  
Total Power Draw appx. 5amps (operating) 9amps (start-up)

#### Dimensions:

Dimension A (WIDTH) -	8'-5"
Dimension A1 (TOP WIDTH) -	8'-0 7/8"
Dimension B (DEPTH) -	3'-9"
Dimension B1 (TOP DEPTH) -	2'-3 7/8"
Dimension C (HEIGHT) -	5'-1 1/2"
Dimension D (TOP HEIGHT) -	3'-0"
Dimension E (CLR. OPEN'G) -	1'-10 1/2"

#### Standard Features:

- Dual HEPA Fan Filter Units w/ Speed Control
- Durable Welded Tube Steel Construction
- Plastic-Laminated Main Work Surface (Shelf Optional)
- White Polyurethane Paint
- 1/2" Side Panels
- Appx.600# with Laminate Top
- 12' Long, 115v Power Cord

#### Surface Options:

- Stainless Steel, Type 304, #4 Finish
- Perforated Stainless Steel, Type 304, #4 Finish
- ESD Laminate Color \_\_\_\_\_
- 180 Degree Rolled Front Edge
- Power Strip  6  8  15
- Clean Room Compatible

#### General Options:


- Cleanroom Shroud
- Casters/Rollers  W/ Brakes
- Cleanroom Casters  W/ Brakes
- ABS Plastic Modular Drawer System
- Metal File Drawers  Metal Utility Drawers
- Lower Shelf  HALF  FULL



**RDM INDUSTRIAL PRODUCTS, INC.**  
1652 Watson Court  
Milpitas, California 95035  
Ph:(408) 945-8400 Fax(408) 945-8433

## Sample Products

### LFH-102P-8 HORIZONTAL LAMINAR FLOW HEPA WORKSTATION



TOTAL UNIT QTY.: \_\_\_\_\_  
**OVERALL DIMENSIONS: 3'-9"D x 8'-5"W x 5'-10 1/2"H**

Light & HEPA FFU are both listed.  
 Total Power Draw appx. 5amps  
 (operating) 9amps (start-up)

#### Dimensions:

Dimension A (WIDTH) -	8'-5"
Dimension A1 (TOP WIDTH) -	8'-0 7/8"
Dimension B (DEPTH) -	3'-9"
Dimension B1 (TOP DEPTH) -	2'-3 7/8"
Dimension C (HEIGHT) -	5'-10 1/2"
Dimension D (TOP HEIGHT) -	3'-0"
Dimension E (CLR. OPEN'G) -	1'-10 1/2"

#### Standard Features:

- Dual HEPA Fan Filter Units w/ Speed Control
- Durable Welded Tube Steel Construction
- Plastic-Laminated Main Work Surface (Shelf Optional)
- White Polyurethane Paint
- 1/2" Side Panels
- Appx. 600# with Laminate Top
- 12' Long, 115v Power Cord

#### Surface Options:

- Stainless Steel, Type 304, #4 Finish
- Perforated Stainless Steel, Type 304, #4 Finish
- ESD Laminate Color \_\_\_\_\_
- 180 Degree Rolled Front Edge
- Power Strip  6  8  15
- Clean Room Compatible

#### General Options:

- Cleanroom Shroud
- Casters/Rollers  W/ Brakes
- Cleanroom Casters  W/ Brakes
- ABS Plastic Modular Drawer System
- Metal File Drawers  Metal Utility Drawers
- Lower Shelf  HALF  FULL

**RDM** RDM INDUSTRIAL PRODUCTS, INC.  
 1652 Watson Court  
 Milpitas, California 95035  
 Ph:(408) 945-8400 Fax(408) 945-8433



# RDM Industrial Products Inc.

QUALITY • SERVICE • TIMING - phone (877) 777-9130

Model LF-102P Vertical



## Model LF-102P Vertical Clean Air Hood:

### Laminar Flow Hood with Removable Work Surface

- Standard Width Sizes include 36", 48", 60", 72" and 96"
- Customizable (RDM can accommodate custom needs)
- Light & Blower Switch are combined into one so that user knows when unit is in operating mode.
- Fluorescent Light Fixture(s)
- Plugs into standard 115v Outlet
- Elec. Junction Box with 12' cord
- Optional Features (Not Included):
  - Full-Height, Vertical Slat Curtain available
  - ULPA Filter(s) available
  - ESD Work Surface
  - Static Control Ionizing Grid(s) available
  - Power Outlet Strip(s) – 8 & 16 Outlet




Page 7



w w w . r d m - i n d . c o m

## Sample Products

### LF-102P-6 LAMINAR FLOW HEPA WORKSTATION



TOTAL UNIT QTY.: \_\_\_\_\_  
**OVERALL DIMENSIONS: 2'-6"D x 6'-0"W x 7'-1 1/2"H**

Light & HEPA FFU are both listed.  
 Total Power Draw appx. 2.5amps  
 (operating) 4.5amps (start-up)

#### Dimensions:

Dimension A (WIDTH) -	6'-0"
Dimension B (DEPTH) -	2'-6"
Dimension C (HEIGHT) -	7'-1 1/2"
Dimension D (TOP HEIGHT) -	3'-0"
Dimension E (CLR. OPEN'G) -	2'-8 1/2"

#### Standard Features:

- Single HEPA Fan Filter Unit w/ Speed Control
- Durable Welded Tube Steel Construction
- Plastic-Laminated Main Work Surface (Shelf Optional)
- White Polyurethane Paint
- 1/4" Clear Acrylic Rear & Side Panels
- Removable Front Draft Guide (40 mil)
- Appx. 350# with Laminate Top
- 12' Long, 115v Power Cord

#### Surface Options:

- Stainless Steel, Type 304, #4 Finish
- Perforated Stainless Steel, Type 304, #4 Finish
- ESD Laminate Color \_\_\_\_\_
- 180 Degree Rolled Front Edge
- Power Strip  6  8  15
- Clean Room Compatible

#### General Options:

- Cleanroom Shroud
- Specialty Panels:  UV  Opaque
- Casters/Rollers  W/ Brakes
- Cleanroom Casters  W/ Brakes
- ABS Plastic Modular Drawer System
- Metal File Drawers  Metal Utility Drawers
- Full Front Panel
- Lower Shelf  HALF  FULL




**RDM INDUSTRIAL PRODUCTS, INC.**  
 1652 Watson Court  
 Milpitas, California 95035  
 Ph:(408) 945-8400 Fax:(408) 945-8433



## Sample Products

### LF-102P-6 LAMINAR FLOW HEPA WORKSTATION



TOTAL UNIT QTY.: \_\_\_\_\_  
**OVERALL DIMENSIONS: 2'-6"D x 6'-0"W x 7'-4 1/2"H**

Light & HEPA FFU are both listed.  
 Total Power Draw appx. 2.5amps  
 (operating) 4.5amps (start-up)

#### Dimensions:

Dimension A (WIDTH) -	6'-0"
Dimension B (DEPTH) -	2'-6"
Dimension C (HEIGHT) -	7'-4 1/2"
Dimension D (TOP HEIGHT) -	3'-0"
Dimension E (CLR. OPEN'G) -	2'-8 1/2"

#### Standard Features:

- Single HEPA Fan Filter Unit w/ Speed Control
- Durable Welded Tube Steel Construction
- Plastic-Laminated Main Work Surface (Shelf Optional)
- White Polyurethane Paint
- 1/4" Clear Acrylic Rear & Side Panels
- Removable Front Draft Guide (40 mil)
- Appx. 350# with Laminate Top
- 12' Long, 115v Power Cord

#### Surface Options:

- Stainless Steel, Type 304, #4 Finish
- Perforated Stainless Steel, Type 304, #4 Finish
- ESD Laminate Color \_\_\_\_\_
- 180 Degree Rolled Front Edge
- Power Strip  6  8  15
- Clean Room Compatible

#### General Options:


- Cleanroom Shroud
- Specialty Panels:  UV  Opaque
- Casters/Rollers  W/ Brakes
- Cleanroom Casters  W/ Brakes
- ABS Plastic Modular Drawer System
- Metal File Drawers  Metal Utility Drawers
- Full Front Panel
- Lower Shelf  HALF  FULL

**RDM** RDM INDUSTRIAL PRODUCTS, INC.  
 1652 Watson Court  
 Milpitas, California 95035  
 Ph:(408) 945-8400 Fax(408) 945-8433

GENERAL NOTE: This drawing is the property of RDM Industrial Products and any use of this drawing without the written consent of RDM is prohibited. No responsibility is assumed by RDM for deviations from this drawing.

## Sample Products

### LF-102P-6-FC2X3 LAMINAR FLOW HEPA WORKSTATION



TOTAL UNIT QTY.: \_\_\_\_\_  
 OVERALL DIMENSIONS: 2'-6"D x 6'-5"W x 7'-1 1/2"H

Light & HEPA FFU are both listed.  
 Total Power Draw appx. 2.5amps  
 (operating) 4.5amps (start-up)

#### Dimensions:

Dimension A (WIDTH) -	6'-5"
Dimension B (DEPTH) -	2'-6"
Dimension C (HEIGHT) -	7'-1 1/2"
Dimension D (TOP HEIGHT) -	3'-0"
Dimension E (CLR. OPEN'G) -	2'-8 1/2"

#### Standard Features:

- Single HEPA Fan Filter Unit w/ Speed Control
- Durable Welded Tube Steel Construction
- Plastic-Laminated Main Work Surface (Shelf Optional)
- White Polyurethane Paint
- 1/4" Clear Acrylic Rear & Side Panels
- Removable Front Draft Guide (40 mil)
- Appx. 350# with Laminate Top
- 12' Long, 115v Power Cord

#### Surface Options:

- Stainless Steel, Type 304, #4 Finish
- Perforated Stainless Steel, Type 304, #4 Finish
- ESD Laminate Color \_\_\_\_\_
- 180 Degree Rolled Front Edge
- Power Strip  6  8  15
- Clean Room Compatible

#### General Options:


- Cleanroom Shroud
- Specialty Panels:  UV  Opaque
- Casters/Rollers  W/ Brakes
- Cleanroom Casters  W/ Brakes
- ABS Plastic Modular Drawer System
- Metal File Drawers  Metal Utility Drawers
- Full Front Panel
- Lower Shelf  HALF  FULL



**RDM INDUSTRIAL PRODUCTS, INC.**  
 1652 Watson Court  
 Milpitas, California 95035  
 Ph:(408) 945-8400 Fax(408) 945-8433

## Sample Products

### LF-102P-8 LAMINAR FLOW HEPA WORKSTATION



TOTAL UNIT QTY.: \_\_\_\_\_

OVERALL DIMENSIONS: 2'-6"D x 8'-5"W x 7'-2" (No Shroud) / 7'-4 1/2"H (w/ Shroud)

Light & HEPA FFU are both listed.  
Total Power Draw appx. 5amps  
(operating) 9amps (start-up)

#### Dimensions:

Dimension A (WIDTH) -	8'-5"
Dimension B (DEPTH) -	2'-6"
Dimension C1/C2 (HEIGHT) -	7'-2" / 7'-4 1/2"
Dimension D (TOP HEIGHT) -	3'-0"
Dimension E (CLR. OPEN'G) -	2'-8 1/2"

#### Standard Features:

- Double HEPA Fan Filter Units w/ Speed Control
- Durable Welded Tube Steel Construction
- Plastic-Laminated Main Work Surface (Shelf Optional)
- White Polyurethane Paint
- 1/4" Clear Acrylic Rear & Side Panels
- Removable Front Draft Guide (40 mil)
- Appx. 450# with Laminate Top
- 12' Long, 115v Power Cord

#### Surface Options:

- Stainless Steel, Type 304, #4 Finish
- Perforated Stainless Steel, Type 304, #4 Finish
- ESD Laminate Color \_\_\_\_\_
- 180 Degree Rolled Front Edge
- Power Strip  6  8  15
- Clean Room Compatible

#### General Options:

- Cleanroom Shroud
- Specialty Panels:  UV  Opaque
- Casters/Rollers  W/ Brakes
- Cleanroom Casters  W/ Brakes
- ABS Plastic Modular Drawer System
- Metal File Drawers  Metal Utility Drawers
- Full Front Panel
- Lower Shelf  HALF  FULL



**RDM INDUSTRIAL PRODUCTS, INC.**  
1652 Watson Court  
Milpitas, California 95035  
Ph: (408) 945-8400 Fax: (408) 945-8433



# RDM Industrial Products Inc.

QUALITY • SERVICE • TIMING - phone (877) 777-9130



## Model LF-102P-B Bench-Top Unit Vertical Clean Air Hood:


A Great Addition to Your Existing Tables...

- Standard Width Sizes include 36", 48", 60", 72" and 96"
- Customizable (RDM can accommodate custom needs)
- Light & Blower Switch are combined into one so that user knows when unit is in operating mode.
- Fluorescent Light Fixture(s)
- Plugs into standard 115v Outlet
- Elec. Junction Box with 12' cord
- Optional Features (Not Included):
  - Full-Height, Vertical Slat Curtain available
  - ULPA Filter(s) available
  - ESD Work Surface
  - Static Control Ionizing Grid(s) available
  - Power Outlet Strip(s) – 8 & 16 Outlet



## Sample Products

### LF-102P-6 LAMINAR FLOW HEPA WORKSTATION - BENCHTOP MODEL



Light & HEPA FFU are both listed.  
Total Power Draw appx. 2.5amps  
(operating) 4.5amps (start-up)

TOTAL UNIT QTY.: \_\_\_\_  
OVERALL DIMENSIONS: 2'-6"D x 6'-0"W x 4'-4 1/2"H

#### Dimensions:

Dimension A (WIDTH) -	6'-0"
Dimension B (DEPTH) -	2'-6"
Dimension C (HEIGHT) -	4'-4 1/2"
Dimension D (CLR. OPEN'G) -	2'-11 1/2"

#### Standard Features:

- Single HEPA Fan Filter Unit w/ Speed Control
- Durable Welded Tube Steel Construction
- White Polyurethane Paint
- 1/4" Clear Acrylic Rear & Side Panels
- Removable Front Draft Guide (40 mil)
- Appx. 225#
- 12' Long, 115v Power Cord

#### ESD Package Options:

- ESD BASE Static Control Package Includes:
  - Plug-in Jacks at Left & Right Uprights
  - ESD Soft Vinyl Front Curtain
  - Clear View Static Control Side Panels
- ESD Static Control Overhead Air Ionizer Grid


#### General Options:

- Cleanroom Shroud
- Specialty Panels:  UV  Opaque
- Double-Sided Access
- Metal File Drawers  Metal Utility Drawers
- Full Front Panel
- Power Strip  6  8  16
- Clean Room Compatible

**RDM** RDM INDUSTRIAL PRODUCTS, INC.  
1652 Watson Court  
Milpitas, California 95035  
Ph:(408) 945-8400 Fax(408) 945-8433

## Sample Products

### LF-102P-6 LAMINAR FLOW HEPA WORKSTATION - BENCHTOP MODEL



Light & HEPA FFU are both listed.  
Total Power Draw appx. 2.5amps  
(operating) 4.5amps (start-up)

TOTAL UNIT QTY.: \_\_\_\_\_  
OVERALL DIMENSIONS: 2'-6"D x 6'-0"W x 4'-7 1/2"H

#### Dimensions:

Dimension A (WIDTH) -	6'-0"
Dimension B (DEPTH) -	2'-6"
Dimension C (HEIGHT) -	4'-7 1/2"
Dimension D (CLR. OPEN'G) -	2'-11 1/2"

#### Standard Features:

- Single HEPA Fan Filter Unit w/ Speed Control
- Durable Welded Tube Steel Construction
- White Polyurethane Paint
- 1/4" Clear Acrylic Rear & Side Panels
- Removable Front Draft Guide (40 mil)
- Appx. 225#
- 12' Long, 115v Power Cord

#### ESD Package Options:

- ESD BASE Static Control Package Includes:
  - Plug-in Jacks at Left & Right Uprights
  - ESD Soft Vinyl Front Curtain
  - Clear View Static Control Side Panels
- ESD Static Control Overhead Air Ionizer Grid

#### General Options:

- Cleanroom Shroud
- Specialty Panels:  UV  Opaque
- Double-Sided Access
- Metal File Drawers  Metal Utility Drawers
- Full Front Panel
- Power Strip  6  8  16
- Clean Room Compatible



**RDM INDUSTRIAL PRODUCTS, INC.**  
1652 Watson Court  
Milpitas, California 95035  
Ph:(408) 945-8400 Fax(408) 945-8433



# RDM Industrial Products Inc.

QUALITY • SERVICE • TIMING - phone (877) 777-9130




Fig. 1 - Front View




Fig. 2 - Rear View

## Instructions:

1. Open main electrical junction box (Fig. 3) located on top of fan filter unit.
2. Record wiring connections before disconnecting and reconnecting new fan filter unit.
3. Disconnect wiring for power supply, light and switch from main electrical box.
4. Remove the (2) Holding Straps at Each end of the Filter (See Fig.2)
5. Replace old fan filter unit with new unit.
6. Secure new fan filter unit into place by replacing the (2) straps (See Fig. 2).
7. Reconnect wiring.

NOTE: Under normal use, assuming eight(8) hours of operation per day, in a standard office environment, fan filter units could last appx. Five (5) years.




Fig. 3 - Electrical Junction Box





# RDM Industrial Products Inc.

QUALITY • SERVICE • TIMING - phone (877) 777-9130


Fan Filter Replacement

**ENVIRCO**  
Innovators in clean air technology

MAC 10 Original

## MAC 10 Original

Extra-Low Fan Filter Unit



### MAC 10 Original

**Extra low watts, extra low sound, extra low profile.**

As the newest addition to Enviroco's MAC 10 Fan Filter Unit (FFU) family, the MAC 10 Original is the only FFU combining low sound, low watts and a low profile. Measuring only 50 dBA, the MAC 10 Original Standard (2 ft x 4 ft [600 mm x 1,210 mm]) provides one of the lowest sound levels of any FFU in the industry. Running at only 200 watts @90 FPM (.45 m/s [175 Watts for 220V/50 Hz]), the unit uses less energy than traditional FFUs, lowering operating costs. The MAC 10 Original maintains a low profile, measuring less than 13 in. (330 mm).

#### STANDARD FEATURES

- Low Sound, Low Watts, Low Profile and Low Operating Costs
- Three Speed Switch features low, medium and high settings; standard on all 2 ft x 4 ft (600 mm x 1,210 mm) and 2 ft x 3.5 ft (600 mm x 1,057 mm) units
- Solid-State Speed Control standard on 2 ft x 2 ft (600 mm x 600 mm) and 2 ft x 3 ft (600 mm x 905 mm) units
- High Efficiency Particulate Air (HEPA) Filter: 99.99% efficient @ 0.3 micron
- Snap-In Prefilter allows for easy replacement and maintenance
- 0.4 in wg (100 Pa) External Static Capability at 90 ft/min (0.45 m/s) with HEPA Filter
- Walkable Plenum (excluding prefilter)
- Meets IEST recommended practices
- UL listed (120 V, 240 V, 277 V)

#### OPTIONS

- Solid-State Speed Control: available on 2 ft x 4 ft (600 mm x 1,210 mm) and 2 ft x 3.5 ft (600 mm x 1,057 mm) units; allows for a full range of settings
- Room-Side Replacement (RSR) available with gel seal filter element (RSR and RSRE models)
- Ultra Low Penetration Air (ULPA) Filter: 99.9995% efficient @ 0.12 micron
- Airflow Indicator Light allows external monitoring of unit operation
- CE Marked: 230V units available
- A/C Collar (12 in [305 mm] diameter) allows direct connection to the air conditioning supply
- Fluorescent Light provides illumination with minimal airflow disruption
- Ion Bar neutralizes static charges below the filter
- Exterior: available powder coat painted or stainless steel
- Power Cord provides convenient plug-in connection for non-permanent applications (non-UL)
- Challenge Port offers convenient aerosol challenge and filter testing
- Knife Edge permits easy placement in gel track ceiling grid system
- Custom Sizes and Configurations available: perfect for mini-environment applications
- Metric Sizes available










# RDM Industrial Products Inc.

QUALITY • SERVICE • TIMING - phone (877) 777-9130


## MAC 10 Original



## MAC 10 Original

**Dimensions**  
inches (millimeters)

- 2 ft x 2ft 23.63 x 23.63 (600 x 600)
- 2 ft x 3 ft 23.63 x 35.63 (600 x 905)
- 2 ft x 3.5 ft 23.63 x 41.63 (600 x 1,057)
- 2 ft x 4 ft 23.63 x 47.63 (600 x 1,210)



MAC 10 Original HEPA Model, 1/5 Hp Motor	AVERAGE AIRFLOW ft <sup>3</sup> /min (m <sup>3</sup> /hr)	HIGH SPEED 120V, 220V, 277V/ 60 Hz (Airflow/Watts)	MEDIUM SPEED 120V, 220V, 277V/ 60 Hz (Airflow/Watts)	LOW SPEED 120V, 220V, 277V/ 60 Hz (Airflow/Watts)	90 ft/min (Airflow/Watts)	70 ft/min (Airflow/Watts)	WATTS 60 Hz (50 Hz)**	SENSIBLE HEAT (BTU/H) 50 Hz @ 90 ft/min (50 Hz @ 0.45 m/s)	SENSIBLE HEAT (BTU/H) 60 Hz @ Full Load (50 Hz @ Full Load)	UNIT WEIGHT lbs (kg)
--	---	--	--	---	---------------------------	---------------------------	-----------------------	--	---	----------------------

### Standard (STD) Model

2 ft x 4 ft 600mm x 1210 mm	650 (1105)	115/270	100/230	80/190	90/200	70/180	200 (180)	683 (615)	922 (829)	66 (29.94)
2 ft x 3.5 ft 600mm x 1057 mm	560 (952)	130/270	115/230	90/190	90/189	70/164	186 (170)	645 (581)	922 (829)	63 (28.58)
2 ft x 3 ft 600mm x 905 mm	480 (816)	150/270			90/163	70/142	163 (147)	556 (501)	922 (829)	52 (23.59)
2 ft x 2 ft 600mm x 600 mm	360 (612)	230/270			90/120		120 (108)	410 (369)	922 (829)	41 (18.6)

### Room-Side Replaceable (RSR) Model

2 ft x 4 ft 600mm x 1210 mm	650 (1,105)	105/270	90/230	70/190	90/220	70/160	220(198)	751 (676)	922 (829)	71 (32.31)
2 ft x 3.5 ft 600mm x 1057 mm	560 (952)	120/270	105/230	80/190	90/180	70/160	180(162)	614 (553)	922 (829)	68 (30.84)
2 ft x 3 ft 600mm x 905 mm	480 (816)	135/270			90/168	70/150	168(151)	573 (517)	922 (829)	54 (24.49)
2 ft x 2 ft 600mm x 600 mm	380 (612)	145/270			90/130	70/120	130(117)	444 (400)	922 (829)	44 (19.96)

Full-Load Amp 120 V / 60 Hz / 1 Ph = 3.0; 220 V / 50/60 Hz / 1 Ph = 1.5; 277 V / 60 Hz / 1 Ph = 1.2

\*\*All test data is measured at 90 ft/min (0.45 m/s). Test Data for the 2 ft x 4 ft (600 mm x 1,210 mm) and 2 ft x 3.5 ft (600 mm x 1,057 mm) measured at 90 ft/min (0.45 m/s) with optional Speed Control; 2 ft x 3 ft (600 mm x 905 mm) and 2 ft x 2 ft (600 mm x 600 mm) units are standard with Speed Control. Test data taken with an electronic micromanometer with Velgrid attachment; eight readings in unloaded conditions. Measurement with two velgrid readings centered on each half of the 2x4 filter gives readings equal to the eight readings.

**ENVIRCO**  
The Greenleaf Centre  
Balloon Fiesta Parkway NE  
Albuquerque, New Mexico 87113  
USA  
Tel: 505.545.3561  
Fax: 505.545.8875  
Toll-Free: 800.545.6598  
email: info@envirco.com

**Trion, Ltd.**  
The Greenleaf Centre  
D - 50171 Kirtson  
Winchester, Hampshire,  
SO22 8JL, UK  
Phone: +44 (0) 1962 840605  
Fax: +44 (0) 1962 828619  
www.trion.co.uk  
email: info@trion.co.uk

**Trion Deutschland GmbH**  
Boschstr. 69  
D - 50171 Kirtson  
New Castle  
Hampshire, Hampshire,  
UK  
Phone: +49 (0) 2227 922-100  
Fax: +49 (0) 2227 922-104  
email: info@trion-deutschland.de

**Feeders International Air Conditioning**  
A-30, Chertsey Road,  
Kallang Colony  
New Delhi  
India 110048  
Phone: (91-11) 5173 0720  
Fax: 91 11 5173 0743  
email: info@indiafeeders.com

**Feeders International**  
No. 728 Futaba Road  
Rm. 3002-3005  
Raffles City Office Tower  
Singapore, P.O. Box, 200001  
Phone: (96-21) 5340 2222  
Fax: (96-21) 5340 3857/5340 3868  
email: info@feeders.com

**FEEDERS**  
Environmental Filtration Systems  
Environmental Filtration Systems  
New Mexico, USA. Feeders is a registered  
manufacturer of products from Enviroco, Inc.  
Specialization in indoor air conditioning  
and energy efficient filtration. ISO certified  
4200, ISO and 14001. Other patents  
issued and pending in foreign countries.

www.envirco.com • 1.800.545.6598

