


Fume Hood w/ Built-in Blower

# RDM Industrial Products Inc.

1652 Watson Court  
 Milpitas, California 95035  
 phone (408) 945-8400 fax (408) 945-8433



**Model B-Max 3/4/5/600:**  
**Fume Hood with Built-in Blower Motor**  
 Provides easy installation for smaller fume control needs...

- Chain-Driven Laminated Safety Glass Sash
- Light & Blower Switches
- Fluorescent Light Fixture
- (2) Total 115v Outlets
- Elec. Junction Box with 7' cord
- Optional Remote Service Fixtures
- 3'-0" & 4'-0" Widths Standard

●\*\* Larger sizes also available, but not recommended because of increased noise due to larger blower motors that are required.

Model #	Exterior Dimensions			Interior Dimensions			2-Tube Light Fixture	Base Cabinets	Shipping Weight
	Width	Depth	Height	Width	Depth	Collar Size			
B-Max-300	36"	33"	59"	27"	26"	9" dia.	1-24" wide	1-36"	375 lbs
B-Max-400	48"	33"	59"	39"	26"	9" dia.	1-36" wide	1-48"	475 lbs
B-Max-500	60"	33"	59"	51"	26"	9" dia.	1-48" wide	2-30"	650 lbs
B-Max-600	72"	33"	59"	63"	26"	9" dia.	1-48" wide	2-36"	725 lbs

