

W-MAX Walk-In Fume Hood

RDM Industrial Products Inc.

1652 Watson Court
 Milpitas, California 95035
 phone (408) 945-8400 fax (408) 945-8433



Model W-Max Fume Hood:




Fume Hood with Large Sliding Doors

Provides more vertical clearance to accommodate your larger project needs...

- Chemical-Resistant Phenolic Liner
- Removable Sliding Doors
- Light & Blower Switches
- Fluorescent Light Fixture(s)
- (2) Total 115v Outlets
- Elec. Junction Box with 7' cord
- Optional Remote Service Fixtures
- Optional Liners, Including Polypropylene, Stainless Steel, PVC, Fiberglass...

GENERAL NOTES:

- FUME HOOD IS MFG. BY GENE SCIENTIFIC INC.
- FUME HOOD LINER MATERIAL IS A 1/4" WHITE TRESPA.
- FUME HOOD, FLUORESCENT LIGHT FIXTURE COMES COMPLETE WITH BULBS.



MODEL #	EXTERIOR DIMENSIONS			INTERIOR		DUCT SIZE
	WIDTH	DEPTH	HEIGHT	WIDTH	DEPTH	
GE4WIFH	48"	36.5"	95"	36"	24"	10"
GE5WIFH	60"	36.5"	95"	48"	24"	10"
GE6WIFH	72"	36.5"	95"	60"	24"	12"
GE8WIFH	96"	36.5"	95"	72"	24"	TWO 10"
GE10WIFH	120"	36.5"	95"	96"	24"	TWO 12"
CUSTOM SIZES AVAILABLE						

