



# STARSYS™ SYSTEM Product Selection Guide



Starsys™ System – Mobile WorkCenters



SXMWS36EOH3



SXMWD36OH3



SXMWT36EOH2



SXMWTK36EBS2



SXMWSK36EBS2

Metro Starsys™ flexible lab furniture is a mobile-based, highly-adaptive system, that enhances lab productivity and space efficiency by empowering the end user with the ability to easily reconfigure and move units around the lab.

The Starsys Modular System creates flexible space. No matter what you need — open, closed, short, tall, mobile units, stationary elements, preconfigured workcenters — Starsys has the elements for you.

With Starsys, you choose from an array of modular units to fit the space and work specifications. It's that easy. And that flexible.

## Benefits:

- Facilitates fast, efficient product or furniture installation.
- The ability to easily reconfigure lab furniture layout as needs change.
- The ability to easily reconfigure unit design.
- Flexibility to design a storage solution to exact needs.
- Corrosion-proof polymer components for a contamination-free product.
- Durable, easy-to-clean, polymer enclosures that won't chip, peel, dent, or rust — ever.
- Rigid anodized aluminum and coated metal substructure, combined with advanced polymers, provides a robust and stable work environment for even the most sensitive equipment.
- A high level of organization.
- Easy access to supplies.
- Greater storage density in a smaller footprint.



**InterMetro Industries Corporation**  
 North Washington Street, Wilkes-Barre, PA 18705  
 Phone: 570-825-2741 • Fax: 570-825-2852  
 For Product Information Call: 1-800-433-2232

L03-252  
 Printed in U.S.A. 1/04  
 Information and specifications are subject to change without notice. Please confirm at time of order.

Copyright © 2004 InterMetro Industries Corp.

46-01



## Specifications

Starsys™ epoxy resin worksurfaces are generally 1" (25mm) thick. Unless otherwise specified, worksurfaces are monolithic and molded from a modified epoxy resin, with a uniform mixture throughout the top forming a thermosetting material. Back of side splashes are shipped separately and affixed on site to epoxy top. Worksurfaces have a smooth, non-glare finish.

Laminate top material is 1 1/8" (29mm) Chip Core with HPL. Laminate thickness is .04" (1.0mm). Phenolic backer thickness is .02" (0.5mm).

Refer to Starsys chemical resistance report for a complete description of Starsys material specifications. A copy of the report can be found at [www.metro.com](http://www.metro.com) or by contacting your local Metro representative.

## Drawers

### Standard Drawers

Description	Inside Dimensions Height/Width/Length		Weight Capacity		Cat. No.
	(in.)	(mm)	(lbs.)	(kg)	
3" (76mm)	2 <sup>3</sup> / <sub>8</sub> x16 <sup>3</sup> / <sub>8</sub> x15 <sup>7</sup> / <sub>8</sub>	60x416x403	25	11.4	SXS3
6" (152mm)	5 <sup>3</sup> / <sub>8</sub> x16 <sup>3</sup> / <sub>8</sub> x15 <sup>7</sup> / <sub>8</sub>	136x416x403	40	18.2	SXS6
9" (230mm)	8 <sup>3</sup> / <sub>8</sub> x16 <sup>3</sup> / <sub>8</sub> x15 <sup>7</sup> / <sub>8</sub>	213x416x403	40	18.2	SXS9

### Heavy-Duty Drawers

Description	Inside Dimensions Height/Width/Length		Weight Capacity							Cat. No.
	(in.)	(mm)	Mobile WorkCenters with Overhead Cabinets		Stationary Applications		Mobile WorkCenters without Overhead Cabinets			
			(lbs.)	(kg)	(lbs.)	(kg)	(lbs.)	(kg)		
<b>Single-Wide Drawers</b>										
3" (76mm)	2 <sup>1</sup> / <sub>4</sub> x16 <sup>3</sup> / <sub>8</sub> x19 <sup>3</sup> / <sub>16</sub>	57x416x487	150	68	150	68	70	32	SXS3HD	
6" (152mm)	5 <sup>1</sup> / <sub>4</sub> x16 <sup>3</sup> / <sub>8</sub> x19 <sup>3</sup> / <sub>16</sub>	133x416x487	150	68	150	68	70	32	SXS6HD	
9" (230mm)	8 <sup>1</sup> / <sub>4</sub> x16 <sup>3</sup> / <sub>8</sub> x19 <sup>3</sup> / <sub>16</sub>	209x416x487	150	68	150	68	70	32	SXS9HD	
12" (305mm)*	11 <sup>1</sup> / <sub>4</sub> x16 <sup>3</sup> / <sub>8</sub> x19 <sup>3</sup> / <sub>16</sub>	285x416x487	150	68	150	68	70	32	SXS12HD	
<b>Double-Wide Drawers</b>										
3" (76mm)	2 <sup>1</sup> / <sub>8</sub> x36 <sup>1</sup> / <sub>4</sub> x19 <sup>3</sup> / <sub>16</sub>	54x920x487	150	68	150	68	110	50	SXD3HD	
6" (152mm)	5 <sup>1</sup> / <sub>8</sub> x36 <sup>1</sup> / <sub>4</sub> x19 <sup>3</sup> / <sub>16</sub>	130x920x487	150	68	150	68	110	50	SXD6HD	
9" (230mm)	8 <sup>1</sup> / <sub>8</sub> x36 <sup>1</sup> / <sub>4</sub> x19 <sup>3</sup> / <sub>16</sub>	206x920x487	150	68	150	68	110	50	SXD9HD	
12" (305mm)	11 <sup>1</sup> / <sub>8</sub> x36 <sup>1</sup> / <sub>4</sub> x19 <sup>3</sup> / <sub>16</sub>	283x920x487	150	68	150	68	110	50	SXD12HD	

\*Note: Single wide 12" (305mm) file drawers include hanging file accessory.

To specify color for drawer pull, order following catalog number for each drawer:

Description	Cat. No.
Dark Taupe	SXDP-TP
Slate Blue	SXDP-BL
Jade	SXDP-GR
Mauve	SXDP-RD
Code Blue	SXDP-CB
Yellow	SXDP-YL
Black	SXDP-BK
Red	SXDP-RE

Note: 2 drawer pulls are required for double-wide heavy-duty drawers.

### Drawer and Door Label Kit

10 per package. 1" (25mm) H. x 16" (406mm) L.  
Cat. No. **SXLABKIT**



# STARSYS™ SYSTEM

## Overhead Cabinets\*

Cabinets include inner panels and mounting rails.  
Refer to shaded door catalog numbers (pages 4&5) for recommended overhead doors.

Description	Overall Dimensions		Height/Depth/Width		Cat. No.
	(in.)	(mm)	(in.)	(mm)	
<b>Single-Wide Cabinets</b>					
24" (610mm)	26 <sup>3</sup> / <sub>8</sub> x13 <sup>1</sup> / <sub>2</sub> x21 <sup>1</sup> / <sub>2</sub>	670x343x546	24 <sup>5</sup> / <sub>8</sub> x12x17	625x341x432	SXSOH24P
27" (686mm)	29 <sup>3</sup> / <sub>8</sub> x13 <sup>1</sup> / <sub>2</sub> x21 <sup>1</sup> / <sub>2</sub>	746x343x546	27 <sup>5</sup> / <sub>8</sub> x12x17	702x341x432	SXSOH27P
30" (760mm)	32 <sup>3</sup> / <sub>8</sub> x13 <sup>1</sup> / <sub>2</sub> x21 <sup>1</sup> / <sub>2</sub>	822x343x546	30 <sup>5</sup> / <sub>8</sub> x12x17	778x341x432	SXSOH30P
33" (838mm)	35 <sup>3</sup> / <sub>8</sub> x13 <sup>1</sup> / <sub>2</sub> x21 <sup>1</sup> / <sub>2</sub>	898x343x546	33 <sup>5</sup> / <sub>8</sub> x12x17	854x341x432	SXSOH33P
36" (915mm)	38 <sup>3</sup> / <sub>8</sub> x13 <sup>1</sup> / <sub>2</sub> x21 <sup>1</sup> / <sub>2</sub>	973x343x546	36 <sup>5</sup> / <sub>8</sub> x12x17	930x341x432	SXSOH36P
48" (1220mm)	50 <sup>3</sup> / <sub>8</sub> x13 <sup>1</sup> / <sub>2</sub> x21 <sup>1</sup> / <sub>2</sub>	1279x343x546	48 <sup>5</sup> / <sub>8</sub> x12x17	1235x341x432	SXSOH48P

<b>Double-Wide Cabinets</b>					
24" (610mm)	26 <sup>3</sup> / <sub>8</sub> x13 <sup>1</sup> / <sub>2</sub> x41 <sup>3</sup> / <sub>8</sub>	670x343x1051	24 <sup>5</sup> / <sub>8</sub> x10x36 <sup>3</sup> / <sub>4</sub>	625x254x933	SXDOH24P
27" (686mm)	29 <sup>3</sup> / <sub>8</sub> x13 <sup>1</sup> / <sub>2</sub> x41 <sup>3</sup> / <sub>8</sub>	746x343x1051	27 <sup>5</sup> / <sub>8</sub> x10x36 <sup>3</sup> / <sub>4</sub>	702x254x933	SXDOH27P
30" (760mm)	32 <sup>3</sup> / <sub>8</sub> x13 <sup>1</sup> / <sub>2</sub> x41 <sup>3</sup> / <sub>8</sub>	822x343x1051	30 <sup>5</sup> / <sub>8</sub> x10x36 <sup>3</sup> / <sub>4</sub>	778x254x933	SXDOH30P
33" (838mm)	35 <sup>3</sup> / <sub>8</sub> x13 <sup>1</sup> / <sub>2</sub> x41 <sup>3</sup> / <sub>8</sub>	898x343x1051	33 <sup>5</sup> / <sub>8</sub> x10x36 <sup>3</sup> / <sub>4</sub>	854x254x933	SXDOH33P
36" (915mm)	38 <sup>3</sup> / <sub>8</sub> x13 <sup>1</sup> / <sub>2</sub> x41 <sup>3</sup> / <sub>8</sub>	973x343x1051	36 <sup>5</sup> / <sub>8</sub> x10x36 <sup>3</sup> / <sub>4</sub>	930x254x933	SXDOH36P
48" (1220mm)	50 <sup>3</sup> / <sub>8</sub> x13 <sup>1</sup> / <sub>2</sub> x41 <sup>3</sup> / <sub>8</sub>	1279x343x1051	48 <sup>5</sup> / <sub>8</sub> x10x36 <sup>3</sup> / <sub>4</sub>	1235x254x933	SXDOH48P

<b>Shelves for Single-Wide Cabinets</b>					
Wire Shelf	1 <sup>3</sup> / <sub>8</sub> x11x17 <sup>15</sup> / <sub>16</sub>	35x279x456			SXSOHWE
Poly Shelf	1 <sup>1</sup> / <sub>2</sub> x12x18	38x305x457			SXSOHPS

Load rating for single-wide shelves, both wire and poly, is 60 lbs. (27kg).

<b>Shelves for Double-Wide Cabinets</b>					
Wire Shelf	1 <sup>3</sup> / <sub>8</sub> x11x36 <sup>11</sup> / <sub>16</sub>	35x279x932			SXDOHWE
Poly Shelf	1 <sup>1</sup> / <sub>2</sub> x12x37 <sup>5</sup> / <sub>8</sub>	38x305x956			SXDOHPS

Load rating for double-wide shelves, both wire and poly, is 120 lbs. (54kg).

\*Triple-wide overhead cabinets are only offered as part of Mobile WorkCenter offering. Contact your Metro representative for more information.

## Overhead Accessories

Description	Height/Depth/Width		Cat. No.
	(in.)	(mm)	
Overhead Light for Single-Wide Cabinets	1.50x5.82x19.75	38x148x502	SXSOHL
Overhead Light for Double-Wide Cabinets	1.50x5.82x39.55	38x148x1005	SXDOHL
Overhead Light for Triple-Wide Cabinets	1.50x5.82x59.35	38x148x1507	SXTOHL
Overhead Light (No Valance)	1.32x4.63x14.25	34x118x362	SXOHLT
Single-Wide Rails (pair)	2.75x.312x21.375	69x8x542	SXOHWB-22
Double-Wide Rails (pair)	2.75x.312x41.188	69x8x1046	SXOHWB-42
Overhead Filler Strip	2.875x.125x49.75	73x3x1264	SXOHFLR-49
Overhead Filler Strip	2.875x.125x63.00	73x3x1600	SXOHFLR-63

## Drawer and Door Label Kit

10 per package. 1" (25mm) H. x 16" (406mm) L. Cat. No. **SXLABKIT**



**Single-Wide Overhead Cabinet with Solid Door**



**Single-Wide Overhead Cabinet with Wire Shelves (2) and Locking, Clear Door**



**Double-Wide Overhead Cabinet with Polymer Shelves and Locking, Clear Door**









# Preconfigured Mobile WorkCenters with Black Epoxy Tops

Description	Cat. No.	Overall Dimensions Width/Length/Height	Cat. No. With Overhead	Overall Dimensions Width/Length/Height
<b>Single-Wide Units</b>				
Single-Wide Heavy-Duty 12" (305mm) (2)	SXMWS33EBS1	28.6x23.6x33.5		
Single-Wide with Standard Drawers	SXMWS33EBS2	28.6x23.6x33.5		
Single-Wide with Standard Drawers	SXMWS36EBS2	28.6x23.6x36.5	SXMWS36EOH2	28.6x23.6x81.5
Single-Wide with Locking Drawers and Poly Shelves	SXMWS36EBS3	28.6x23.6x36.5	SXMWS36EOH3	28.6x23.6x81.5
Single-Wide with Locking Drawers and Wire Shelves	SXMWS36EBS4	28.6x23.6x36.5	SXMWS36EOH4	28.6x23.6x81.5
Single-Wide with Standard Drawers	SXMWS39EBS2	28.6x23.6x39.5	SXMWS39EOH2	28.6x23.6x81.5
Single-Wide with Locking Drawers and Poly Shelves	SXMWS39EBS3	28.6x23.6x39.5	SXMWS39EOH3	28.6x23.6x81.5
Single-Wide with Locking Drawers and Wire Shelves	SXMWS39EBS4	28.6x23.6x39.5	SXMWS39EOH4	28.6x23.6x81.5
<b>Double-Wide Units</b>				
Double-Wide Heavy-Duty 12" (305mm) (2)	SXMWD33EBS1	28.6x43.4x33.5		
Double-Wide with Standard Drawers	SXMWD33EBS2	28.6x43.4x33.5		
Double-Wide with Standard Drawers	SXMWD36EBS2	28.6x43.4x36.5	SXMWD36EOH2	28.6x43.4x81.5
Double-Wide with DHD Drawers	SXMWD36EBS5	28.6x43.4x36.5	SXMWD36EOH5	28.6x43.4x81.5
Double-Wide with Locking Drawers and Poly Shelves	SXMWD36EBS3	28.6x43.4x36.5	SXMWD36EOH3	28.6x43.4x81.5
Double-Wide with Locking Drawers and Wire Shelves	SXMWD36EBS4	28.6x43.4x36.5	SXMWD36EOH4	28.6x43.4x81.5
Double-Wide with Standard Drawers	SXMWD39EBS2	28.6x43.4x39.5	SXMWD39EOH2	28.6x43.4x81.5
Double-Wide with DHD Drawers	SXMWD39EBS5	28.6x43.4x39.5	SXMWD39EOH5	28.6x43.4x81.5
Double-Wide with Locking Drawers and Poly Shelves	SXMWD39EBS3	28.6x43.4x39.5	SXMWD39EOH3	28.6x43.4x81.5
Double-Wide with Locking Drawers and Wire Shelves	SXMWD39EBS4	28.6x43.4x39.5	SXMWD39EOH4	28.6x43.4x81.5
<b>Triple-Wide Units</b>				
Triple-Wide with Standard Drawers	SXMWT36EBS2	28.6x63.2x36.5	SXMWT36EOH2	28.6x63.2x81.5
Triple-Wide with S&D HD Drawers	SXMWT36EBS6	28.6x63.2x36.5	SXMWT36EOH6	28.6x63.2x81.5
Triple-Wide with Locking Drawers and Poly Shelves	SXMWT36EBS3	28.6x63.2x36.5	SXMWT36EOH3	28.6x63.2x81.5
Triple-Wide with Locking Drawers and Wire Shelves	SXMWT36EBS4	28.6x63.2x36.5	SXMWT36EOH4	28.6x63.2x81.5
Triple-Wide with Standard Drawers	SXMWT39EBS2	28.6x63.2x39.5	SXMWT39EOH2	28.6x63.2x81.5
Triple-Wide with S&D HD Drawers	SXMWT39EBS6	28.6x63.2x39.5	SXMWT39EOH6	28.6x63.2x81.5
Triple-Wide with Locking Drawers and Poly Shelves	SXMWT39EBS3	28.6x63.2x39.5	SXMWT39EOH3	28.6x63.2x81.5
Triple-Wide with Locking Drawers and Wire Shelves	SXMWT39EBS4	28.6x63.2x39.5	SXMWT39EOH4	28.6x63.2x81.5
Triple Kneewell with Standard Drawers	SXMWTK36EBS2	28.6x63.2x36.5	SXMWTK36EOH2	28.6x63.2x81.5
Triple Kneewell with Locking Drawers and Poly Shelves	SXMWTK36EBS3	28.6x63.2x36.5	SXMWTK36EOH3	28.6x63.2x81.5
Triple Kneewell with Locking Drawers and Wire Shelves	SXMWTK36EBS4	28.6x63.2x36.5	SXMWTK36EOH4	28.6x63.2x81.5
Triple Kneewell with Standard Drawers	SXMWTK39EBS2	28.6x63.2x39.5	SXMWTK39EOH2	28.6x63.2x81.5
Triple Kneewell with Locking Drawers and Poly Shelves	SXMWTK39EBS3	28.6x63.2x39.5	SXMWTK39EOH3	28.6x63.2x81.5
Triple Kneewell with Locking Drawers and Wire Shelves	SXMWTK39EBS4	28.6x63.2x39.5	SXMWTK39EOH4	28.6x63.2x81.5
Single Kneewell with HD Drawers	SXMWSK33EBS1	28.6x73x33.5		
Single Kneewell with Standard Drawers	SXMWSK36EBS2	28.6x73x36.5		
Single Kneewell with Locking Drawers and Poly Shelves	SXMWSK36EBS3	28.6x73x36.5		
Single Kneewell with Locking Drawers and Wire Shelves	SXMWSK36EBS4	28.6x73x36.5		
Single Kneewell with Standard Drawers	SXMWSK39EBS2	28.6x73x39.5		
Single Kneewell with Locking Drawers and Poly Shelves	SXMWSK39EBS3	28.6x73x39.5		
Single Kneewell with Locking Drawers and Wire Shelves	SXMWSK39EBS4	28.6x73x39.5		

**Note:** Above units include taupe door and/or drawer pulls.

**Note:** Mobile units include 5" (127mm) poly plate casters (2 toe brake, 2 swivel). Refer to next page for other caster options. Overhead units include a 27" (686mm) overhead cabinet with locking clear doors and two poly shelves. If additional poly shelves are required, please use Catalog Number shown on page 3. All tops are black epoxy resin. Colored laminate and other top materials are available on request. Contact your Metro representative for more information.

# Preconfigured Mobile WorkCenter Units with Laminate Tops

Description	Cat. No.	Overall Dimensions Width/Length/Height
Single-Wide, Locking with Standard Drawers and Overhead	SXSENT1	28.6x23.6x81.5
Double-Wide, Locking with Standard Drawers, Poly Shelf and Overhead	SXDENT2	28.6x43.4x81.5
Double-Wide, Locking with Standard Drawers and Overhead	SXD40MW10	28.6x43.4x81.5
Double-Wide, Locking with Overhead and Wire Shelf	SXD40MW12	28.6x43.4x81.5
Double-Wide, Locking with Standard Drawers, Poly Shelf and Overhead	SXD40MW13	28.6x43.4x81.5
Triple-Wide, with Standard Drawers, Wire Shelf and Overhead	SXT38MW14	28.6x63.2x79.6
Single-Wide Kneewell with Standard and Heavy-Duty Drawers	SXK32MW5	28.6x73x31.9
Single-Wide Kneewell, Locking with Standard Drawers	SXK32MW6	28.6x73x31.9
Single-Wide Kneewell with Standard Drawers and Wire Shelf	SXK41MW7	28.6x73x40.3
Single-Wide Kneewell with Poly Shelf	SXK41MW8	28.6x73x40.3

**Note:** Above units, where appropriate, include a slate blue countertop along with slate blue door and/or drawer pulls.

**Note:** Laminate top material is 1/8" (29mm) Chip Core with HPL. Laminate thickness is .04" (1.0mm). Phenolic backer thickness is .02" (0.5mm).

# STARSYS™ SYSTEM



A sampling of available preconfigured units is shown below. For a complete preconfigured unit overview visit [www.metro.com](http://www.metro.com) or contact your local representative.



SXMWS36EBS2



SXMWD36EBS2



SXMWT36EBS6



SXMWT36EBS3



SXMWS36EBS3



SXMWSK36EBS2



SXMWT36EOH2



SXMWTK36EOH2



SXMWD36EOH3

**Caster Options:** 5" (127mm) polyurethane plate casters come standard with preconfigured units (shown above) and make units easy to maneuver in tight spaces. Other caster options are listed below.

**3" (76mm) Total-Lock Caster:** Utilize 3" (76mm) total-lock plate casters for limited mobile applications or when additional storage area is required.



Cat. No. B3P-TL

**5" (127mm) All-Polymer Total-Lock Caster:** Recommended for corrosive environments. All-polymer total-lock plate caster is economical alternative to stainless casters.



Cat. No. B5PC-TL

**Stabilizer/Leveling Caster:** Recommended for applications where benchtop equipment demands unit stability and/or a level worksurface. Engage caster foot for stability, release for mobility.



Cat. No. B3N-AS

**Total # of Units:** \_\_\_\_\_  
 (fill out rest of WorkSheet based on one unit)  
**Pull Color Options**  
 \*specify color at the unit level  
 Taupe \_\_\_ Slate Blue \_\_\_ Jade \_\_\_ Mauve \_\_\_  
 Red \_\_\_ Black \_\_\_ Code Blue \_\_\_ Yellow \_\_\_

**Unit Lock Options**

Mechanical Lock \_\_\_\_\_  
 For Multiple Locks indicate keyed:  
 Alike \_\_\_ Differently \_\_\_  
 Key # if needed: \_\_\_\_\_

**Polymer Interior Options**

**Polymer Drawers & Accessories**  
 3 Dividers \_\_\_\_\_  
 Short \_\_\_ Long \_\_\_ Kit \_\_\_\_\_  
 6"/9" Dividers \_\_\_\_\_  
 Short \_\_\_ Long \_\_\_ Kit \_\_\_\_\_  
 Label Kits (10/pk) \_\_\_\_\_

**Cassettes & Accessories**

2 Level (9" total height) \_\_\_\_\_  
 3 Level (12" total height) \_\_\_\_\_  
 4 Level (15" total height) \_\_\_\_\_  
 5 Level (18" total height) \_\_\_\_\_  
 3" Bin (6 per level) \_\_\_\_\_  
 4.5" Bin (4 per level) \_\_\_\_\_  
 6" Bin (3 per level) \_\_\_\_\_  
 Cassette Label Kit \_\_\_\_\_  
 (18 bins/size/color) \_\_\_\_\_  
 Bin Divider Kit (up to 4 bins) \_\_\_\_\_

**Poly Shelf Assemblies & Accessories**

Single wide Poly Shelf \_\_\_\_\_  
 Poly Shelf Divider Kit \_\_\_\_\_  
 Poly Shelf w/Drawer Face \_\_\_\_\_

**qwikSLOT Interior Options**

**HD Drawer Accessories**  
 3" Div Kit \_\_\_\_\_ Single \_\_\_ Double \_\_\_  
 6"/9" Div Kit \_\_\_\_\_ Single \_\_\_ Double \_\_\_  
 Label Kits (10 pk) \_\_\_\_\_  
**Shelf Assemblies & Accessories**  
 Drawer-Top Module: \_\_\_\_\_ Single \_\_\_ Double \_\_\_  
 qwikSLOT Shelves: \_\_\_\_\_ Single \_\_\_ Double \_\_\_  
 Full Ext. Shelves: \_\_\_\_\_ Single \_\_\_ Double \_\_\_  
 Wire Shelf Dividers: (21" x 3") \_\_\_\_\_  
 Wire Shelf Ledges: (15 1/2" x 4 1/2") \_\_\_\_\_  
 Label Holders: \_\_\_\_\_ Single \_\_\_ Double \_\_\_

**Mark Starsys Drawers:**

"LK" for locking cover (polymer drawers only)  
 "HD" for heavy-duty drawers  
 "PL" for tamper evident passive locks

**Caster Options**

3" Leveling \_\_\_\_\_  
 3" Total Lock \_\_\_\_\_  
 5" Standard \_\_\_\_\_  
 5" Poly \_\_\_\_\_

**Starsys Door Options**

**Lower Doors**  
 Solid Clear \_\_\_\_\_  
 90° 270° Locking \_\_\_\_\_  
 Left Right \_\_\_\_\_  
**Upper Doors**  
 Solid Clear \_\_\_\_\_  
 90° 270° Locking \_\_\_\_\_  
 Left Right \_\_\_\_\_

**Side Pods & Accessories**

Stationary Side Storage Unit Left Right \_\_\_\_\_  
 Tambour Door \_\_\_\_\_  
 Tilt out/Removable Bins \_\_\_\_\_  
 Full Width Shelf \_\_\_\_\_  
 Top Tray w/dividers \_\_\_\_\_  
 Chart Holder \_\_\_\_\_  
 Gas Tank Cylinder Holder \_\_\_\_\_  
 Sharps Container Holder \_\_\_\_\_  
 Cup Holder \_\_\_\_\_  
 Waste Basket & Holder \_\_\_\_\_

**Tambour Door Overhead**

Overhead Unit \_\_\_\_\_ Quantity \_\_\_\_\_  
 Upper Shelf \_\_\_\_\_  
 MB30230B (12 per carton) \_\_\_\_\_  
 MB30235B (6 per carton) \_\_\_\_\_  
 3" Bin (6 per shelf) \_\_\_\_\_  
 4.5" Bin (4 per shelf) \_\_\_\_\_  
 6" Bin (3 per shelf) \_\_\_\_\_  
 Tilt Bin \_\_\_\_\_  
 Lower Shelf \_\_\_\_\_  
 Under Cabinet Light \_\_\_\_\_

**Hinged Door Overhead**

(half depth Starsys cabinets)  
 Height: 24" 27" 30" 33" 36" 48"  
 Single-wide Double-wide \_\_\_\_\_  
 Poly Shelves Qty. \_\_\_\_\_  
 Wire Shelves Qty. \_\_\_\_\_  
 Poly Shelf Div Kit Qty. \_\_\_\_\_  
 Single Double \_\_\_\_\_  
 Under Cabinet Light: \_\_\_\_\_  
 Slanted Tops: (for wall mount only)  
 Single Double Triple \_\_\_\_\_

**Mobile WorkCenter Countertops**

Black Epoxy Laminated \_\_\_\_\_  
 Laminate Colors \_\_\_\_\_  
 Slate Blue Mauve Jade \_\_\_\_\_  
 Custom Color Spec: \_\_\_\_\_

**Counter Details**

Back splash \_\_\_\_\_ Side splash \_\_\_\_\_  
 OH Supports \_\_\_\_\_ Island Top \_\_\_\_\_

**Kneewell Options**

Structural Support Options \_\_\_\_\_  
 Pencil Drawer 30" 42"  
 Keyboard Tray 30" 42"



Worksheet for:  
**Mobile WorkCenters, Mobile WorkCenters with Overhead Cabinets**



Worksheet \_\_\_\_\_ of \_\_\_\_\_ Date: \_\_\_\_\_ Account #: \_\_\_\_\_

Account Name: \_\_\_\_\_ Project Name: \_\_\_\_\_

Send to: \_\_\_\_\_

Fax #: \_\_\_\_\_ E-mail: \_\_\_\_\_

Rep \_\_\_\_\_ Territory #: \_\_\_\_\_



\*Note: Cassettes and 3" Drawers cannot be placed in the bottom 3" of unit.