



# RDM INDUSTRIAL PRODUCTS, INC.


Laboratory & Industrial Furnishing Solutions...




QUALITY • SERVICE • TIMING

# TABLE OF CONTENTS

Specifications-----	1
Base Cabinets - Standing/Counter Height-----	2 - 8
Base Cabinets - Sink/Fume Hood Standing Height-----	9
Base Cabinets - Special Purpose-----	10
Acid Storage Cabinets-----	11
Flammable Storage Cabinets-----	11
Base Cabinets - Sitting/Desk Height-----	12 - 14
Corner Base Cabinets - Sitting/Desk Height-----	15
Corner Sink Cabinets - Sitting/Desk Height-----	15
Wall Case Units-----	16 - 17
Corner Wall Case Units-----	17
Floor Case Units-----	18 - 19
Fume Hoods-----	20 - 21
Aprons and Drawer Housings-----	22
Utility Tables-----	23
Accessories-----	24 - 25



# BASE CABINETS *Standing/Counter Height*



## Additional Specifications

All surfaces shall be painted or zinc plated. The paint is a chemical resistant baked on epoxy powder. Specific test data is available upon request.

Gauges of base cabinet components are as follows:

Aprons	18 ga.	Legs - 2" sq. tube	18 ga.
Back Panels	20 ga.	Shelves	18 ga.
Bottom Panels	18 ga.	Side Panels	18 ga.
Drawer/Door Outer Panel	20 ga.	Table Frames	18 ga.
Door Inner Panel	20 ga.	K & V Shelf Clips	
Drawer Bodies	20 ga.		

NOTE: Cabinet Heights are 29" (Desk-Height) & 35" (Standard Height). Levelers extend 1-1/2".

See ADA-Height Metal laboratory Cabinets catalogue for 32.5" high cabinets.

Subject to change.



# BASE CABINETS *Standing/Counter Height*

Standing/Counter Height 35"High (34.75"w/out Leveler) x 22"Deep. Units are available in various widths as shown. Drawer units feature full extension slides. A single shelf is standard in the units shown. Additional shelves are available (see


**Part Numbering System**

Sample:	RB35H36W-048	1 = Base unit	4 = #Doors
	11222333-456	2 = Height	5 = #Drawers
		3 = Width	6 = #Half Drawers


- \*"F" = File Drawer
- \*"S" = Sink Unit
- \*"P" = Panel
- \*"R" = Right hinge facing cab.
- \*"L" = Left hinge facing cab.




**RB35H15W-000  
RB35H18W-000  
RB35H24W-000**




**RB35H15W-100(L)  
RB35H18W-100(L)  
RB35H24W-100(L)**




**RB35H15W-100(R)  
RB35H18W-100(R)  
RB35H24W-100(R)**




**RB35W15H-110(L)  
RB35W18H-110(L)  
RB35W24H-110(L)**



**RB35H15W-110(R)  
RB35H18W-110(R)  
RB35H24W-110(R)**




**RB35H15W-102(L)  
RB35H18W-102(L)  
RB35H24W-102(L)**




**RB35H15W-102(R)  
RB35H18W-102(R)  
RB35H24W-102(R)**

**NOTE:** When ordering units with configurations that are non-standard, drawings must accompany order.


**BASE CABINETS** *Standing/Counter Height*




**RB35H15W-F20  
RB35H18W-F20  
RB35H24W-F20**




**RB35H15W-040  
RB35H18W-040  
RB35H24W-040**




**RB35H15W-032  
RB35H18W-032  
RB35H24W-032**




**RB35H15W-024  
RB35H18W-024  
RB35H24W-024**




**RB35H15W-008  
RB35H18W-008  
RB35H24W-008**



**RB35H15W-F12  
RB35H18W-F12  
RB35H24W-F12**




**RB35H15W-FFO  
RB35H18W-FFO  
RB35H24W-FFO**




**RB35H30W-000  
RB35H36W-000  
RB35H42W-000  
RB35H48W-000**


**BASE CABINETS** *Standing/Counter Height*




**RB35H30W-200  
RB35H36W-200  
RB35H42W-200  
RB35H48W-200**




**RB35H30W-110(L)SP  
RB35H36W-110(L)SP  
RB35H42W-110(L)SP  
RB35H48W-110(L)SP**




**RB35H30W-110(R)SP  
RB35H36W-110(R)SP  
RB35H42W-110(R)SP  
RB35H48W-110(R)SP**  
NOTE: 1/2 CAB+1" Exposed




**RB35H30W-210  
RB35H36W-210  
RB35H42W-210  
RB35H48W-210**




**RB35H30W-220  
RB35H36W-220  
RB35H42W-220  
RB35H48W-220**



**RB35H30W-212A  
RB35H36W-212A  
RB35H42W-212A  
RB35H48W-212A**




**RB35H30W-212B  
RB35H36W-212B  
RB35H42W-212B  
RB35H48W-212B**




**RB35H30W-204  
RB35H36W-204  
RB35H42W-204  
RB35H48W-204**


**BASE CABINETS** *Standing/Counter Height*




RB35H30W-140B(L)  
 RB35H36W-140B(L)  
 RB35H42W-140B(L)  
 RB35H48W-140B(L)




RB35H30W-140B(R)  
 RB35H36W-140B(R)  
 RB35H42W-140B(R)  
 RB35H48W-140B(R)




RB35H30W-150(L)  
 RB35H36W-150(L)  
 RB35H42W-150(L)  
 RB35H48W-150(L)




RB35H30W-150(R)  
 RB35H36W-150(R)  
 RB35H42W-150(R)  
 RB35H48W-150(R)




RB35H30W-142(L)  
 RB35H36W-142(L)  
 RB35H42W-142(L)  
 RB35H48W-142(L)



RB35H30W-142(R)  
 RB35H36W-142(R)  
 RB35H42W-142(R)  
 RB35H48W-142(R)




RB35H30W-132(L)  
 RB35H36W-132(L)  
 RB35H42W-132(L)  
 RB35H48W-132(L)




RB35H30W-132(R)  
 RB35H36W-132(R)  
 RB35H42W-132(R)  
 RB35H48W-132(R)


**BASE CABINETS** *Standing/Counter Height*




RB35H30W-100(L)SP  
 RB35H36W-100(L)SP  
 RB35H42W-100(L)SP  
 RB35H48W-100(L)SP




RB35H30W-100(R)SP  
 RB35H36W-100(R)SP  
 RB35H42W-100(R)SP  
 RB35H48W-100(R)SP




RB35H30W-140A(L)  
 RB35H36W-140A(L)  
 RB35H42W-140A(L)  
 RB35H48W-140A(L)




RB35H30W-140A(R)  
 RB35H36W-140A(R)  
 RB35H42W-140A(R)  
 RB35H48W-140A(R)




RB35H30W-134(L)  
 RB35H36W-134(L)  
 RB35H42W-134(L)  
 RB35H48W-134(L)



RB35H30W-134(R)  
 RB35H36W-134(R)  
 RB35H42W-134(R)  
 RB35H48W-134(R)




RB35H30W-048A  
 RB35H36W-048A  
 RB35H42W-048A  
 RB35H48W-048A




RB35H30W-048B  
 RB35H36W-048B  
 RB35H42W-048B  
 RB35H48W-048B

**NOTE:** When ordering units with configurations that are non-standard, drawings must accompany order.


# BASE CABINETS *Standing/Counter Height*




**RB35H30W-040**  
**RB35H36W-040**  
**RB35H42W-040**  
**RB35H48W-040**  
**NOTE: Drawer I.D.:**  
**Width=CABWidth-4.875**  
**Depth=19.5"D**




**RB35H30W-050**  
**RB35H36W-050**  
**RB35H42W-050**  
**RB35H48W-050**




**RB35H30W-060**  
**RB35H36W-060**  
**RB35H42W-060**  
**RB35H48W-060**




**RB35H30W-080**  
**RB35H36W-080**  
**RB35H42W-080**  
**RB35H48W-080**




**RB35H30W-024**  
**RB35H36W-024**  
**RB35H42W-024**  
**RB35H48W-024**



**RB35H30W-032**  
**RB35H36W-032**  
**RB35H42W-032**  
**RB35H48W-032**




**RB35H30W-044A**  
**RB35H36W-044A**  
**RB35H42W-044A**  
**RB35H48W-044A**




**RB35H30W-044B**  
**RB35H36W-044B**  
**RB35H42W-044B**  
**RB35H48W-044B**


**BASE CABINETS** *Standing/Counter Height*




**RB35H30W-048C  
RB35H36W-048C  
RB35H42W-048C  
RB35H48W-048C**




**RB35H30W-064  
RB35H36W-064  
RB35H42W-064  
RB35H48W-064**




**RB35H30W-070A  
RB35H36W-070A  
RB35H42W-070A  
RB35H48W-070A**



**RB35H30W-070B  
RB35H36W-070B  
RB35H42W-070B  
RB35H48W-070B**



**RB35H30W-008  
RB35H36W-008  
RB35H42W-008  
RB35H48W-008**



**RB35H30W-052  
RB35H36W-052  
RB35H42W-052  
RB35H48W-052**

**NOTE:** When ordering units with configurations that are non-standard, drawings must accompany order.

# SINK AND FUME HOOD BASE CABINETS *Stand Up Height*

ALSO AVAILABLE IN SITTING/DESK HEIGHT



**Sink Base**  
RB35H18W-10S(R)  
RB35H24W-10S(R)



**Sink Base**  
RB35H18W-10S(L)  
RB35H24W-10S(L)



**Sink Base**  
RB35H18W-10S(R)  
Louvered  
RB35H24W-10S(R)  
Louvered




**Sink Base**  
RB35H18W-10S(L)  
Louvered  
RB35H24W-10S(L)  
Louvered




**Hood Base**  
RB35H18W-10P(R)  
RB35H24W-10P(R)




**Hood Base**  
RB35H18W-10P(L)  
RB35H24W-10P(L)



**Sink Base**  
RB35H30W-20S  
RB35H36W-20S  
RB35H42W-20S  
RB35H48W-20S  
RB35H52W-20S  
RB35H58W-20S



**Sink Island Base**  
RB35H48W-20SI  
RB35H52W-20SI  
RB35H58W-20SI




**Sink Base**  
RB35H30W-20S Louvered  
RB35H36W-20S Louvered  
RB35H42W-20S Louvered  
RB35H48W-20S Louvered  
RB35H52W-20S Louvered  
RB35H58W-20S Louvered



**Hood Base**  
RB35H30W-20P  
RB35H36W-20P  
RB35H42W-20P  
RB35H48W-20P  
RB35H52W-20P  
RB35H58W-20P




**Hood Island Base**  
RB35H52W-20PI  
RB35H58W-20PI



**Sink Support Assembly**  
SS-36  
SS-C14 SS-48  
SS-18 SS-42  
SS-24 SS-52

**BASE CABINETS** *Corner Unit Bases – w/ Removable Backs*



Typical Dimensions

**RC35H14W-100(L)**



**RC35H14W-100(R)**  
*(SHOWN WITH  
OPTIONAL LAZY  
SUSAN INSTALLED)*



**RC35H14W-10S(L)**



**RC35H14W-10S(R)**




**RC35H14W-10S(L)**  
Louvered




**RC35H14W-10S(R)**  
Louvered

**MOBILE BASES**



**RMB29H15W-030**  
**RMB29H18W-030**  
**RMB29H24W-030**  
**RMB29H30W-030**




**RMB35H15W-040**  
**RMB35H18W-040**  
**RMB35H24W-040**  
**RMB35H30W-040**

**NOTE: 3" Casters (4" O.H.) can be added to Any Cabinet Style in lieu of Toe Kick Base. Be sure to add back panel if creating mobile modular pedestal style cabinet.**

# ACID STORAGE CABINET

Metal-Arc Acid Storage Cabinets are reinforced steel with a hi-tech performance polymer alloy coating. This reduces the possibilities of corrosion. Note: All Acid Cabinets are 18" deep.




- |               |                |
|---------------|----------------|
| RASB29H18-100 | RASB35H18W-100 |
| RASB29H24-100 | RASB35H24W-100 |
| RASB29H30-200 | RASB35H30W-200 |
| RASB29H36-200 | RASB35H36W-200 |
| RASB29H42-200 | RASB35H42W-200 |
| RASB29H48-200 | RASB35H48W-200 |
| RASB29H52-200 | RASB35H52W-200 |
| RASB29H58-200 | RASB35H58W-200 |

} 1 Door Units

Add Suffix "VK" if Vent Kit is Required  
Add Suffix "SCC" if Self-Closing Doors Required

# FLAMMABLE STORAGE CABINET

Note: All Flammable Cabinets are 18" deep.




- |                |                |
|----------------|----------------|
| RFSB29H24W-100 | RFSB35H24W-100 |
| RFSB29H30W-200 | RFSB35H30W-200 |
| RFSB29H36W-200 | RFSB35H36W-200 |
| RFSB29H42W-200 | RFSB35H42W-200 |
| RFSB29H48W-200 | RFSB35H48W-200 |
| RFSB29H52W-200 | RFSB35H52W-200 |
| RFSB29H58W-200 | RFSB35H58W-200 |

} 1 Door Units

# STANDARD SHELF

Standard base cabinet shelf with brackets. Available for 29" and 35" units.




- RBS15W
- RBS18W
- RBS21W
- RBS24W
- RBS30W
- RBS36W
- RBS42W
- RBS48W

# PULL-OUT SHELF

Optional 2" high pull-out shelf for base cabinets. Shelf is full extension with ball bearing slides. 19.5" deep. \*Note: Cabinet door must be open 90 deg. to allow for extension.

8 max. 35"H cab.  
6 max. 29"H cab.



- RBPS15W
- RBPS18W
- RBPS21W
- RBPS24W
- RBPS30W
- RBPS36W
- RBPS42W
- RBPS48W



NOTE: When ordering units with configurations that are non-standard, drawings must accompany order.

# BASE CABINETS *Sitting / Desk Height*


Sitting/Desk Height 29" High X 22" Deep. Units are available in various widths as shown. Drawer units feature full extension slides. A single shelf is standard in the units where shown. Additional shelves are available (see page 11). Specify left or right hand doors. Hand is indicated by side of door upon which hinges are located.

### Part Numbering System

- |         |              |              |                  |
|---------|--------------|--------------|------------------|
| Sample: | RB29H30W-044 | 1= Base Unit | 4= #Doors        |
|         | 11222333-456 | 2= Height    | 5= #Drawers      |
|         |              | 3= Width     | 6= #Half Drawers |
|         |              |              | 7= Features      |
- \* "F" = File Drawer
  - \* "S" = Sink Unit
  - \* "P" = Panel over door(s)
  - \* "R" = Right hinge facing cab.
  - \* "L" = Left hinge facing cab.
  - \* "B" = Bin Drawer




- RB29H15W-000
- RB29H18W-000
- RB29H24W-000
- RB29H30W-000
- RB29H36W-000
- RB29H42W-000
- RB29H48W-000




### End View


- RB29H15W-100(L)
- RB29H18W-100(L)
- RB29H24W-100(L)




- RB29H15W-100(R)
- RB29H18W-100(R)
- RB29H24W-100(R)




- RB29H15W-110(L)
- RB29H18W-110(L)
- RB29H24W-110(L)




- RB29H15W-110(R)
- RB29H18W-110(R)
- RB29H24W-110(R)



- RB29H15W-102(L)
- RB29H18W-102(L)
- RB29H24W-102(L)




- RB29H15W-102(R)
- RB29H18W-102(R)
- RB29H24W-102(R)




- RB29H15W-F10
- RB29H18W-F10
- RB29H24W-F10

**NOTE:** When ordering units with configurations that are non-standard, drawings must accompany order.


# BASE CABINETS *Sitting / Desk Height*




RB29H15W-F02  
RB29H18W-F02  
RB29H24W-F02




RB29H30W-040  
RB29H36W-040  
RB29H42W-040  
RB29H48W-040




RB29H30W-050  
RB29H36W-050  
RB29H42W-050  
RB29H48W-050




RB29H30W-044  
RB29H36W-044  
RB29H42W-044  
RB29H48W-044




RB29H30W-006  
RB29H36W-006  
RB29H42W-006  
RB29H48W-006



RB29H30W-060  
RB29H36W-060  
RB29H42W-060  
RB29H48W-060




RB29H30W-200  
RB29H36W-200  
RB29H42W-200  
RB29H48W-200




RB29H30W-210  
RB29H36W-210  
RB29H42W-210  
RB29H48W-210


**BASE CABINETS** *Sitting / Desk Height*




**RB29H15W-030  
RB29H18W-030  
RB29H24W-030**




**RB29H15W-B11  
RB29H18W-B11  
RB29H24W-B11**



**RB29H15W-022  
RB29H18W-022  
RB29H24W-022**



**RB29H15W-006  
RB29H18W-006  
RB29H24W-006**



**RB29H30W-100(RorL)SP  
RB29H36W-100(RorL)SP  
RB29H42W-100(RorL)SP  
RB29H48W-100(RorL)SP**



**RB29H30W-220  
RB29H36W-220  
RB29H42W-220  
RB29H48W-220**




**RB29H30W-110(R)SP  
RB29H36W-110(R)SP  
RB29H42W-110(R)SP  
RB29H48W-110(R)SP**




**RB29H30W-110(L)SP  
RB29H36W-110(L)SP  
RB29H42W-110(L)SP  
RB29H48W-110(L)SP**

**CORNER BASE CABINETS** *Sitting / Desk Height*



Typical Dimensions


**RC29H14W-100(L)**




**RC29H14W-100(R)**

*(SHOWN WITH OPTIONAL LAZY SUSAN INSTALLED)*


**CORNER SINK BASE CABINETS** *Sitting / Desk Height*




**RC29H14W-10S(L)**



**RC29H14W-10S(R)**




**RC29H14W-10S(L)**  
Louvered



**RC29H14W-10S(R)**  
Louvered

**OPTIONAL LAZY SUSAN** *(For corner base cabinets)*



**RLS-29H**  
**RLS-35H**

# WALL CASE UNITS


Wall cases are available in sizes ranging from 24", 31", 36" and 48" high and 18" wide to 48" wide. All units are 13" deep (11.75" deep I.D. w/ 10.5" D Shelves) as standard. 24" high units have one shelf. 31" high to 36" high units have two shelves. 48" high units have three shelves. Additional shelves are available. On single hinge door units, specify (R) or (L) hinge. 48" high units have 3 hinges per door.

## Part Numbering System

Sample: RW30H36W-2GS 1= Wall Case 4= #Doors  
 11222333-4567 2= Height 5= Material of doors\*  
 3= Width 6= Type of doors\* \*


- \* Steel, Glass or no doors (open)
- \*\* Hinged, Sliding or Sliding/framed
- \* "R" = Right hinge facing cabinet
- \* "L" = Left hinge facing cabinet

7= Depth Note: 13" dp. standard,  
 16" deep cabinets available -  
 indicate by final suffix of 16,  
 (Example: RW24H18W-00016)




### OPEN SHELVING UNITS:

RW24H18W-000	RW30H18W-000	RW36H18W-000	RW48H18W-000
RW24H24W-000	RW30H24W-000	RW36H24W-000	RW48H24W-000
RW24H30W-000	RW30H30W-000	RW36H30W-000	RW48H30W-000
RW24H36W-000	RW30H36W-000	RW36H36W-000	RW48H36W-000
RW24H42W-000	RW30H42W-000	RW36H42W-000	RW48H42W-000
RW24H48W-000	RW30H48W-000	RW36H48W-000	RW48H48W-000




### STEEL DOOR HINGED LEFT:

RW24H18W-1SH(L)	RW30H18W-1SH(L)
RW36H18W-1SH(L)	RW48H18W-1SH(L)
RW24H24W-1SH(L)	RW30H24W-1SH(L)
RW36H24W-1SH(L)	RW48H24W-1SH(L)




### STEEL DOOR HINGED RIGHT:

RW24H18W-1SH(R)	RW30H18W-1SH(R)
RW36H18W-1SH(R)	RW48H18W-1SH(R)
RW24H24W-1SH(R)	RW30H24W-1SH(R)
RW36H24W-1SH(R)	RW48H24W-1SH(R)




### GLASS W/ FRAME HINGED LEFT:

RW24H18W-1GH(L)	RW30H18W-1GH(L)
RW36H18W-1GH(L)	RW48H18W-1GH(L)
RW24H24W-1GH(L)	RW30H24W-1GH(L)
RW36H24W-1GH(L)	RW48H24W-1GH(L)




### GLASS W/ FRAME HINGED RIGHT:

RW24H18W-1GH(R)	RW30H18W-1GH(R)
RW36H18W-1GH(R)	RW48H18W-1GH(R)
RW24H24W-1GH(R)	RW30H24W-1GH(R)
RW36H24W-1GH(R)	RW48H24W-1GH(R)



### STEEL DOORS - HINGED:


RW24H30W-2SH	RW36H30W-2SH
RW24H36W-2SH	RW36H36W-2SH
RW24H42W-2SH	RW36H42W-2SH
RW24H48W-2SH	RW36H48W-2SH
RW30H30W-2SH	RW48H30W-2SH
RW30H36W-2SH	RW48H36W-2SH
RW30H42W-2SH	RW48H42W-2SH
RW30H48W-2SH	RW48H48W-2SH



### STEEL DOORS - SLIDING:


RW24H30W-2SS	RW36H30W-2SS
RW24H36W-2SS	RW36H36W-2SS
RW24H42W-2SS	RW36H42W-2SS
RW24H48W-2SS	RW36H48W-2SS
RW30H30W-2SS	RW48H30W-2SS
RW30H36W-2SS	RW48H36W-2SS
RW30H42W-2SS	RW48H42W-2SS
RW30H48W-2SS	RW48H48W-2SS

# WALL CASE UNITS




**GLASS W/ FRAME – HINGED:**

- RW24H30W-2GH W36H30W-2GH
- RW24H36W-2GH W36H36W-2GH
- RW24H42W-2GH W36H42W-2GH
- RW24H48W-2GH W36H48W-2GH
- RW30H30W-2GH W48H30W-2GH
- RW30H36W-2GH W48H36W-2GH
- RW30H42W-2GH W48H42W-2GH
- RW30H48W-2GH W48H48W-2GH




**GLASS DOORS – SLIDING:**

- RW24H30W-2GS W36H30W-2GS
- RW24H36W-2GS W36H36W-2GS
- RW24H42W-2GS W36H42W-2GS
- RW24H48W-2GS W36H48W-2GS
- RW30H30W-2GS W48H30W-2GS
- RW30H36W-2GS W48H36W-2GS
- RW30H42W-2GS W48H42W-2GS
- RW30H48W-2GS W48H48W-2GS




**GLASS W/ FRAME – SLIDING:**

- RW24H30W-2GSF R W36H30W-2GSF
- RW24H36W-2GSF R W36H36W-2GSF
- RW24H42W-2GSF R W36H42W-2GSF
- RW24H48W-2GSF R W36H48W-2GSF
- RW30H30W-2GSF R W48H30W-2GSF
- RW30H36W-2GSF R W48H36W-2GSF
- RW30H42W-2GSF R W48H42W-2GSF
- RW30H48W-2GSF R W48H48W-2GSF



**PASS-THRU (OPEN ON TWO SIDES)**


- RWDS30H30W-00013 open-2-s
- RWDS30H36W-00013 open-2-s
- RWDS30H48W-00013 open-2-s



**PASS-THRU (OPEN ON TWO SIDES)**


- RWDS30H36W-4GS13 open-2-s
- RWDS30H48W-4GS13 open-2-s

# CORNER WALL CASE UNITS




**OPEN SHELVING:**

- RWC24H18W-000
- RWC30H18W-000
- RWC36H18W-000
- RWC30H18W-00016
- RWC36H18W-00016




**STEEL DOOR – HINGED:**

- RWC24H18W-1SH(RorL)
- RWC30H18W-1SH(RorL)
- RWC36H18W-1SH(RorL)
- RWC30H18W-1SH(RorL)16
- RWC36H18W-1SH(RorL)16



**GLASS DOOR W/ FRAME:**

- RWC24H18W-1GH(L)
- RWC30H18W-1GH(L)
- RWC36H18W-1GH(L)
- RWC30H18W-1GH(L)16
- RWC36H18W-1GH(L)16




**GLASS DOOR W/ FRAME:**

- RWC24H18W-1GH(R)
- RWC30H18W-1GH(R)
- RWC36H18W-1GH(R)
- RWC30H18W-1GH(R)16
- RWC36H18W-1GH(R)16

**NOTE:** When ordering units with configurations that are non-standard, drawings must accompany order.  
**NOTE:** 16" deep corner wall case dimension = 28-11/16"




# FLOOR CASE UNITS




**GLASS W/ FRAME – HINGED LEFT:**

- RF84H18W-1GH16(L)
- RF84H18W-1GH21(L)
- RF84H24W-1GH16(L)
- RF84H24W-1GH21(L)




**GLASS W/ FRAME – HINGED RIGHT:**

- RF84H18W-1GH16(R)
- RF84H18W-1GH21(R)
- RF84H24W-1GH16(R)
- RF84H24W-1GH21(R)




**STEEL DOORS – HINGED:**

- RF84H30W-2SH16
- RF84H30W-2SH21
- RF84H36W-2SH16
- RF84H36W-2SH21
- RF84H42W-2SH16
- RF84H42W-2SH21
- RF84H48W-2SH16
- RF84H48W-2SH21




**STEEL DOORS – SLIDING:**

- RF84H30W-2SS16
- RF84H30W-2SS21
- RF84H36W-2SS16
- RF84H36W-2SS21
- RF84H42W-2SS16
- RF84H42W-2SS21
- RF84H48W-2SS16
- RF84H48W-2SS21



**GLASS W/ FRAME – SLIDING:**

- RF84H30W-2GSF16
- RF84H30W-2GSF21
- RF84H36W-2GSF16
- RF84H36W-2GSF21
- RF84H42W-2GSF16
- RF84H42W-2GSF21
- RF84H48W-2GSF16
- RF84H48W-2GSF21



**GLASS W/ FRAME – HINGED:**

- RF84H30W-2GH16
- RF84H30W-2GH21
- RF84H36W-2GH16
- RF84H36W-2GH21
- RF84H42W-2GH16
- RF84H42W-2GH21
- RF84H48W-2GH16
- RF84H48W-2GH21

**NOTE:** When ordering units with configurations that are non-standard, drawings must accompany order.

## FUME HOODS

### QUALITY


The name RDM has long been synonymous with quality and longevity. Our newly redesigned fume hood line is a proud example of this reputation. We are backed with over 50 years of experience, combined with state of the art engineering, manufacturing and testing facilities to provide you with a top quality fume hood solution unsurpassed in quality and performance.

### DESIGN

The design of our fume hood line is founded on three basic principals: safety, flexibility and aesthetics.

### SAFETY

Safety should be the number one concern when choosing a fume hood. That is why our hoods are tested and certified in accordance with the Ashrae method of testing performance of laboratory fume hoods.



### FLEXIBILITY

Construction of our fume hood begins with a heavy duty internal framework. This universal framework is the foundation for a variety of standard and optional features.

These features include. . .

- Chemical resistant interiors of resin-chem, PVC or stainless steel
- PVC or stainless sash guides
- Vertical, horizontal or combination sashes in a variety of materials
- Easy access for service fixture installation
- Manual or remote baffle adjustment
- Fluorescent light fixture
- Optional factory installation of service fixtures

The flexibility of our design will allow you to select standard options that will meet the requirements of your application.

### AESTHETICS

Our quality, pride and craftsmanship will insure that your RDM Fume Hood will be a handsome center piece in your laboratory for years to come.



### CATALOG

For a complete fume hood catalog contact your local RDM representative or contact us direct at the following address:

**RDM Industrial Products, Inc.**  
**2020 Traders Ridge Drive**  
**Conroe, Texas 77301**  
**PH: (936) 777-8700**  
**FX: (936) 777-8707**



# FUME HOODS

**BYPASS & VARIABLE  
AIR VOLUME (VAV)  
FUME HOODS**



**PERCHLORIC ACID AND  
RADIO-ISOTOPE FUME HOODS**


**ADA FUME HOODS**




**WALK-IN AND SPECIAL  
FUME HOODS**

# APRONS AND DRAWER HOUSINGS


Part Numbering System: RA04H60W22D-000 1= Apron 4= Depth  
 11222333444-555 2= Height 5= # Drawers  
 3= Width




- |                 |                 |
|-----------------|-----------------|
| RA04H24W22D-000 | RA04H24W28D-000 |
| RA04H30W22D-000 | RA04H30W28D-000 |
| RA04H36W22D-000 | RA04H36W28D-000 |
| RA04H42W22D-000 | RA04H42W28D-000 |
| RA04H48W22D-000 | RA04H48W28D-000 |
| RA04H54W22D-000 | RA04H54W28D-000 |
| RA04H60W22D-000 | RA04H60W28D-000 |
| RA04H66W22D-000 | RA04H66W28D-000 |
| RA04H72W22D-000 | RA04H72W28D-000 |



- |                 |                 |
|-----------------|-----------------|
| RA04H24W22D-001 | RA04H24W28D-001 |
| RA04H30W22D-001 | RA04H30W28D-001 |
| RA04H36W22D-001 | RA04H36W28D-001 |
| RA04H42W22D-001 | RA04H42W28D-001 |




- |                 |                 |
|-----------------|-----------------|
| RA04H42W22D-002 | RA04H42W28D-002 |
| RA04H48W22D-002 | RA04H48W28D-002 |
| RA04H54W22D-002 | RA04H54W28D-002 |
| RA04H60W22D-002 | RA04H60W28D-002 |
| RA04H66W22D-002 | RA04H66W28D-002 |
| RA04H72W22D-002 | RA04H72W28D-002 |



- |                 |                 |
|-----------------|-----------------|
| RA04H72W22D-003 | RA04H72W28D-003 |
|-----------------|-----------------|

Note: Apron heights are 4-1/2" high.

# KNEE SPACE RAILS




- |            |            |
|------------|------------|
| RKSR04H18W | RKSR04H52W |
| RKSR04H24W | RKSR04H54W |
| RKSR04H30W | RKSR04H60W |
| RKSR04H36W | RKSR04H66W |
| RKSR04H42W | RKSR04H72W |
| RKSR04H48W | RKSR04H84W |
|            | RKSR04H96W |

NOTE: When ordering units with configurations that are non-standard, drawings must accompany order.

# UTILITY TABLES AND COMPONENTS

## 29" HIGH TABLES

## 35" HIGH TABLES



**22" DEEP**

RUT29H24W22D-000  
RUT29H30W22D-000  
RUT29H36W22D-000  
RUT29H42W22D-000  
RUT29H48W22D-000  
RUT29H60W22D-000  
RUT29H72W22D-000

**28" DEEP**


RUT29H24W28D-000  
RUT29H30W28D-000  
RUT29H36W28D-000  
RUT29H42W28D-000  
RUT29H48W28D-000  
RUT29H60W28D-000  
RUT29H72W28D-000

**22" DEEP**

RUT35H24W22D-000  
RUT35H30W22D-000  
RUT35H36W22D-000  
RUT35H42W22D-000  
RUT35H48W22D-000  
RUT35H60W22D-000  
RUT35H72W22D-000

**28" DEEP**

RUT35H24W28D-000  
RUT35H30W28D-000  
RUT35H36W28D-000  
RUT35H42W28D-000  
RUT35H48W28D-000  
RUT35H60W28D-000  
RUT35H72W28D-000



**22" DEEP**

RUT29H24W22D-001  
RUT29H30W22D-001  
RUT29H36W22D-001  
RUT29H42W22D-001

**28" DEEP**


RUT29H24W28D-001  
RUT29H30W28D-001  
RUT29H36W28D-001  
RUT29H42W28D-001

**22" DEEP**

RUT35H24W22D-001  
RUT35H30W22D-001  
RUT35H36W22D-001  
RUT35H42W22D-001

**28" DEEP**

RUT35H24W28D-001  
RUT35H30W28D-001  
RUT35H36W28D-001  
RUT35H42W28D-001



**22" DEEP**

RUT29H48W22D-002  
RUT29H60W22D-002  
RUT29H72W22D-002

**28" DEEP**


RUT29H48W28D-002  
RUT29H60W28D-002  
RUT29H72W28D-002

**22" DEEP**

RUT35H48W22D-002  
RUT35H60W22D-002  
RUT35H72W22D-002

**28" DEEP**

RUT35H48W28D-002  
RUT35H60W28D-002  
RUT35H72W28D-002



**22" DEEP**

RUT29H72W22D-003

**28" DEEP**

RUT29H72W28D-003

**22" DEEP**

RUT35H72W22D-003

**28" DEEP**


RUT35H72W28D-003

**NOTE:** UTILITY TABLES DO NOT INCLUDE COUNTER TOPS. 2" SQ. LEGS, BLACK VINYL BOOTS.

**NOTE:** When ordering units with configurations that are non-standard, **drawings must accompany order.**

# ACCESSORIES

### TABLE LEGS




- RL29H-22D
- RL35H-22D
- RL29H-28D
- RL35H-28D

### SINGLE LEG




- RSL-29H
- RSL-35H

### SINGLE LEG WITH PLATE




- RSL-29HP
- RSL-35HP

### END PEDESTAL




- REP29H02W-22D
- REP35H02W-22D
- REP29H02W-28D
- REP35H02W-28D

### WALL CASE (UPPER CAB) FRONT FILLER



- Specify Width (W)
- RWCFF24H-W
  - RWCFF30H-W
  - RWCFF36H-W
  - RWCFF48H-W

### WALL CASE CORNER FILLER




Specify Depth (D)  
13" or 16"

- RWCCF24H-D
- RWCCF30H-D
- RWCCF36H-D
- RWCCF48H-D

Angle is 1"

### REAGENT SHELF UPRIGHT




- RRSU18H12D2W
- (18"Hx12"Dx2"W)

Note: Special heights & options available upon request. Drawings and dimensions must accompany order.

NOTE: Use with RSF Fascia Panels for Lateral Top Support (Not Shown) - Mounts between Uprights.

### WALL - FLOOR CASE SHELF



### SLOPE TOP WITH END CAPS

#### PART NUMBERS

#### PART NUMBERS


#### WALL

- RWS18W13D
- RWS24W13D
- RWS30W13D
- RWS36W13D
- RWS42W13D
- RWS48W13D
- RWS18W16D
- RWS24W16D
- RWS30W16D
- RWS36W16D
- RWS42W16D
- RWS48W16D

#### FLOOR

- RFS18W16D
- RFS24W16D
- RFS30W16D
- RFS36W16D
- RFS42W16D
- RFS48W16D
- RFS18W21D
- RFS24W21D
- RFS30W21D
- RFS36W21D
- RFS42W21D
- RFS48W21D


- RST13D-18W
- RST13D-24W
- RST13D-30W
- RST13D-36W
- RST13D-42W
- RST13D-48W
- RST16D-18W
- RST16D-24W
- RST16D-30W
- RST16D-36W
- RST16D-42W
- RST16D-48W
- RST21D-18W
- RST21D-24W
- RST21D-30W
- RST21D-36W
- RST21D-42W
- RST21D-48W



NOTE: When ordering units with configurations that are non-standard, drawings must accompany order.

# ACCESSORIES

## KNEE SPACE PANELS



For use in sit down situation with apron

### 24" HIGH

RKSP24H24W  
RKSP24H30W  
RKSP24H36W  
RKSP24H42W  
RKSP24H48W  
RKSP24H54W  
RKSP24H60W  
RKSP24H66W  
RKSP24H72W

For use in standing situation with apron, or sit down without apron

### 29" HIGH

RKSP29H24W  
RKSP29H30W  
RKSP29H36W  
RKSP29H42W  
RKSP29H48W  
RKSP29H54W  
RKSP29H60W  
RKSP29H66W  
RKSP29H72W

For use in standing situation without apron

### 35" HIGH

RKSP35H24W  
RKSP35H30W  
RKSP35H36W  
RKSP35H42W  
RKSP35H48W  
RKSP35H54W  
RKSP35H60W  
RKSP35H66W  
RKSP35H72W


# FILL PANELS

All base cabinet fillers are 35" high. Fillers for 29" high (low base) cabinets will need to be trimmed by the installer. Other sizes available upon request.




### C-FILLER

RCF35-8W  
RCF35-14W



### L-FILLER

RLF35-1W  
RLF35-7W  
RLF35-11W



### C-FRONT FILLER

RFFC35-8W  
RFFC35-14W



### CORNER FILLER

RICF35H-1W  
RICF35H-2W  
RICF35H-4W



### L-FRONT FILLER

RFFL35H-8W  
RFFL35H-14W



### FLOOR CASE FILLER

RFFL84H-2W  
RFFL84H-8W




### TOE SPACE FILLER

RTSF-00 (WIDTH)

\* Last number(s) in filler designation indicates width of filler assembly.



Represented By:



## **RDM** LABORATORY FURNITURE SOLUTIONS

2020 Traders Ridge Drive

Conroe, Texas 77301

Telephone: (877) 777-9130

FAX: (936) 777-8707

Web page: [www.rdm-ind.com](http://www.rdm-ind.com)

E-mail: [lab-sales@rdm-ind.com](mailto:lab-sales@rdm-ind.com)

© RDM Industrial Products, Inc.

2025 Printed in the U.S.A.

**QUALITY • SERVICE • TIMING**