



RDM INDUSTRIAL PRODUCTS, INC.


Laboratory & Industrial Furnishing Solutions...




QUALITY • SERVICE • TIMING

TABLE OF CONTENTS

Specifications-----	1
Base Cabinets - Standing/Counter Height-----	2 - 8
Base Cabinets - Sink/Fume Hood Standing Height-----	9
Base Cabinets - Special Purpose-----	10
Acid Storage Cabinets-----	11
Flammable Storage Cabinets-----	11
Base Cabinets - Sitting/Desk Height-----	12 - 14
Corner Base Cabinets - Sitting/Desk Height-----	15
Corner Sink Cabinets - Sitting/Desk Height-----	15
Wall Case Units-----	16 - 17
Corner Wall Case Units-----	17
Floor Case Units-----	18 - 19
Fume Hoods-----	20 - 21
Aprons and Drawer Housings-----	22
Utility Tables-----	23
Accessories-----	24 - 25



BASE CABINETS *Standing/Counter Height*



Additional Specifications

All surfaces shall be painted or zinc plated. The paint is a chemical resistant baked on epoxy powder. Specific test data is available upon request.

Gauges of base cabinet components are as follows:

Aprons	18 ga.	Legs - 2" sq. tube	18 ga.
Back Panels	20 ga.	Shelves	18 ga.
Bottom Panels	18 ga.	Side Panels	18 ga.
Drawer/Door Outer Panel	20 ga.	Table Frames	18 ga.
Door Inner Panel	20 ga.	K & V Shelf Clips	
Drawer Bodies	20 ga.		

NOTE: Cabinet Heights are 29" (Desk-Height) & 35" (Standard Height). Levelers extend 1-1/2".

See ADA-Height Metal laboratory Cabinets catalogue for 32.5" high cabinets.

Subject to change.

BASE CABINETS *Standing/Counter Height*

Standing/Counter Height 35"High (34.75"w/out Leveler) x 22"Deep. Units are available in various widths as shown. Drawer units feature full extension slides. A single shelf is standard in the units shown. Additional shelves are available (see

Part Numbering System

Sample:

RB35H36W-048
11222333-456

1 = Base unit
2 = Height
3 = Width

4 = #Doors
5 = #Drawers
6 = #Half Drawers


*"F" = File Drawer

*"S" = Sink Unit


*"P" = Panel

*"R" = Right hinge facing cab.


*"L" = Left hinge facing cab.




**RB35H15W-000
RB35H18W-000
RB35H24W-000**




**RB35H15W-100(L)
RB35H18W-100(L)
RB35H24W-100(L)**




**RB35W15H-110(L)
RB35W18H-110(L)
RB35W24H-110(L)**




**RB35H15W-102(L)
RB35H18W-102(L)
RB35H24W-102(L)**




**RB35H15W-100(R)
RB35H18W-100(R)
RB35H24W-100(R)**



**RB35H15W-110(R)
RB35H18W-110(R)
RB35H24W-110(R)**




**RB35H15W-102(R)
RB35H18W-102(R)
RB35H24W-102(R)**




NOTE: When ordering units with configurations that are non-standard, drawings must accompany order.


BASE CABINETS *Standing/Counter Height*




**RB35H15W-F20
RB35H18W-F20
RB35H24W-F20**




**RB35H15W-040
RB35H18W-040
RB35H24W-040**




**RB35H15W-032
RB35H18W-032
RB35H24W-032**




**RB35H15W-024
RB35H18W-024
RB35H24W-024**




**RB35H15W-008
RB35H18W-008
RB35H24W-008**



**RB35H15W-F12
RB35H18W-F12
RB35H24W-F12**




**RB35H15W-FFO
RB35H18W-FFO
RB35H24W-FFO**



**RB35H30W-000
RB35H36W-000
RB35H42W-000
RB35H48W-000**


BASE CABINETS *Standing/Counter Height*




RB35H30W-200
RB35H36W-200
RB35H42W-200
RB35H48W-200




RB35H30W-110(L)SP
RB35H36W-110(L)SP
RB35H42W-110(L)SP
RB35H48W-110(L)SP




RB35H30W-110(R)SP
RB35H36W-110(R)SP
RB35H42W-110(R)SP
RB35H48W-110(R)SP
NOTE: ½ CAB+1" Exposed




RB35H30W-210
RB35H36W-210
RB35H42W-210
RB35H48W-210




RB35H30W-220
RB35H36W-220
RB35H42W-220
RB35H48W-220



RB35H30W-212A
RB35H36W-212A
RB35H42W-212A
RB35H48W-212A




RB35H30W-212B
RB35H36W-212B
RB35H42W-212B
RB35H48W-212B




RB35H30W-204
RB35H36W-204
RB35H42W-204
RB35H48W-204


BASE CABINETS *Standing/Counter Height*




RB35H30W-140B(L)
RB35H36W-140B(L)
RB35H42W-140B(L)
RB35H48W-140B(L)




RB35H30W-140B(R)
RB35H36W-140B(R)
RB35H42W-140B(R)
RB35H48W-140B(R)




RB35H30W-150(L)
RB35H36W-150(L)
RB35H42W-150(L)
RB35H48W-150(L)




RB35H30W-150(R)
RB35H36W-150(R)
RB35H42W-150(R)
RB35H48W-150(R)




RB35H30W-142(L)
RB35H36W-142(L)
RB35H42W-142(L)
RB35H48W-142(L)



RB35H30W-142(R)
RB35H36W-142(R)
RB35H42W-142(R)
RB35H48W-142(R)




RB35H30W-132(L)
RB35H36W-132(L)
RB35H42W-132(L)
RB35H48W-132(L)




RB35H30W-132(R)
RB35H36W-132(R)
RB35H42W-132(R)
RB35H48W-132(R)


BASE CABINETS *Standing/Counter Height*




RB35H30W-100(L)SP
RB35H36W-100(L)SP
RB35H42W-100(L)SP
RB35H48W-100(L)SP




RB35H30W-100(R)SP
RB35H36W-100(R)SP
RB35H42W-100(R)SP
RB35H48W-100(R)SP




RB35H30W-140A(L)
RB35H36W-140A(L)
RB35H42W-140A(L)
RB35H48W-140A(L)




RB35H30W-140A(R)
RB35H36W-140A(R)
RB35H42W-140A(R)
RB35H48W-140A(R)




RB35H30W-134(L)
RB35H36W-134(L)
RB35H42W-134(L)
RB35H48W-134(L)



RB35H30W-134(R)
RB35H36W-134(R)
RB35H42W-134(R)
RB35H48W-134(R)




RB35H30W-048A
RB35H36W-048A
RB35H42W-048A
RB35H48W-048A




RB35H30W-048B
RB35H36W-048B
RB35H42W-048B
RB35H48W-048B

NOTE: When ordering units with configurations that are non-standard, drawings must accompany order.


BASE CABINETS *Standing/Counter Height*




RB35H30W-040
RB35H36W-040
RB35H42W-040
RB35H48W-040
NOTE: Drawer I.D.:
Width=CABWidth-4.875
Depth=19.5"D




RB35H30W-050
RB35H36W-050
RB35H42W-050
RB35H48W-050




RB35H30W-060
RB35H36W-060
RB35H42W-060
RB35H48W-060




RB35H30W-080
RB35H36W-080
RB35H42W-080
RB35H48W-080




RB35H30W-024
RB35H36W-024
RB35H42W-024
RB35H48W-024



RB35H30W-032
RB35H36W-032
RB35H42W-032
RB35H48W-032




RB35H30W-044A
RB35H36W-044A
RB35H42W-044A
RB35H48W-044A




RB35H30W-044B
RB35H36W-044B
RB35H42W-044B
RB35H48W-044B


BASE CABINETS *Standing/Counter Height*




RB35H30W-048C
RB35H36W-048C
RB35H42W-048C
RB35H48W-048C




RB35H30W-064
RB35H36W-064
RB35H42W-064
RB35H48W-064




RB35H30W-070A
RB35H36W-070A
RB35H42W-070A
RB35H48W-070A



RB35H30W-070B
RB35H36W-070B
RB35H42W-070B
RB35H48W-070B



RB35H30W-008
RB35H36W-008
RB35H42W-008
RB35H48W-008



RB35H30W-052
RB35H36W-052
RB35H42W-052
RB35H48W-052

NOTE: When ordering units with configurations that are non-standard, drawings must accompany order.

SINK AND FUME HOOD BASE CABINETS *Stand Up Height*

ALSO AVAILABLE IN SITTING/DESK HEIGHT



Sink Base
RB35H18W-10S(R)
RB35H24W-10S(R)



Sink Base
RB35H18W-10S(L)
RB35H24W-10S(L)



Sink Base
RB35H18W-10S(R)
Louvered
RB35H24W-10S(R)
Louvered




Sink Base
RB35H18W-10S(L)
Louvered
RB35H24W-10S(L)
Louvered




Hood Base
RB35H18W-10P(R)
RB35H24W-10P(R)




Hood Base
RB35H18W-10P(L)
RB35H24W-10P(L)



Sink Base
RB35H30W-20S
RB35H36W-20S
RB35H42W-20S
RB35H48W-20S
RB35H52W-20S
RB35H58W-20S



Sink Island Base
RB35H48W-20SI
RB35H52W-20SI
RB35H58W-20SI




Sink Base
RB35H30W-20S Louvered
RB35H36W-20S Louvered
RB35H42W-20S Louvered
RB35H48W-20S Louvered
RB35H52W-20S Louvered
RB35H58W-20S Louvered



Hood Base
RB35H30W-20P
RB35H36W-20P
RB35H42W-20P
RB35H48W-20P
RB35H52W-20P
RB35H58W-20P




Hood Island Base
RB35H52W-20PI
RB35H58W-20PI



Sink Support Assembly

SS-C14	SS-36
SS-18	SS-48
SS-24	SS-42
	SS-52

BASE CABINETS *Corner Unit Bases – w/ Removable Backs*



Typical Dimensions

RC35H14W-100(L)



RC35H14W-100(R)
(SHOWN WITH
OPTIONAL LAZY
SUSAN INSTALLED)



RC35H14W-10S(L)



RC35H14W-10S(R)




RC35H14W-10S(L)
Louvered




RC35H14W-10S(R)
Louvered

MOBILE BASES



RMB29H15W-030
RMB29H18W-030
RMB29H24W-030
RMB29H30W-030




RMB35H15W-040
RMB35H18W-040
RMB35H24W-040
RMB35H30W-040

**NOTE: 3" Casters (4" O.H.) can be added to Any Cabinet Style in lieu of Toe Kick Base.
Be sure to add back panel if creating mobile modular pedestal style cabinet.**

ACID STORAGE CABINET

Metal-Arc Acid Storage Cabinets are reinforced steel with a hi-tech performance polymer alloy coating. This reduces the possibilities of corrosion. Note: All Acid Cabinets are 18" deep.



RASB29H18-100	RASB35H18W-100
RASB29H24-100	RASB35H24W-100
RASB29H30-200	RASB35H30W-200
RASB29H36-200	RASB35H36W-200
RASB29H42-200	RASB35H42W-200
RASB29H48-200	RASB35H48W-200
RASB29H52-200	RASB35H52W-200
RASB29H58-200	RASB35H58W-200


} 1 Door Units

Add Suffix "VK" if Vent Kit is Required

Add Suffix "SCC" if Self-Closing Doors Required

FLAMMABLE STORAGE CABINET

Note: All Flammable Cabinets are 18" deep.




RFSB29H24W-100
RFSB29H30W-200
RFSB29H36W-200
RFSB29H42W-200
RFSB29H48W-200
RFSB29H52W-200
RFSB29H58W-200

RFSB35H24W-100
RFSB35H30W-200
RFSB35H36W-200
RFSB35H42W-200
RFSB35H48W-200
RFSB35H52W-200
RFSB35H58W-200

} 1 Door Units

STANDARD SHELF

Standard base cabinet shelf with brackets. Available for 29" and 35" units.




RBS15W
RBS18W
RBS21W
RBS24W
RBS30W
RBS36W
RBS42W
RBS48W

PULL-OUT SHELF

Optional 2" high pull-out shelf for base cabinets. Shelf is full extension with ball bearing slides. 19.5" deep.
 *Note: Cabinet door must be open 90 deg. to allow for extension.

8 max. 35"H cab.
 6 max. 29"H cab.



RBPS15W
RBPS18W
RBPS21W
RBPS24W
RBPS30W
RBPS36W
RBPS42W
RBPS48W

NOTE: When ordering units with configurations that are non-standard, drawings must accompany order.



BASE CABINETS *Sitting / Desk Height*

Sitting/Desk Height 29" High X 22" Deep. Units are available in various widths as shown. Drawer units feature full extension slides. A single shelf is standard in the units where shown. Additional shelves are available (see page 11). Specify left or right hand doors. Hand is indicated by side of door upon which hinges are located.


Part Numbering System

Sample: RB29H30W-044 1= Base Unit 4= #Doors
11222333-456 2= Height 5= #Drawers
3= Width 6= #Half Drawers
7= Features

- * "F" = File Drawer
- * "S" = Sink Unit
- * "P" = Panel over door(s)
- * "R" = Right hinge facing cab.
- * "L" = Left hinge facing cab.
- * "B" = Bin Drawer




RB29H15W-000
RB29H18W-000
RB29H24W-000
RB29H30W-000
RB29H36W-000
RB29H42W-000
RB29H48W-000




End View


RB29H15W-100(L)
RB29H18W-100(L)
RB29H24W-100(L)




RB29H15W-100(R)
RB29H18W-100(R)
RB29H24W-100(R)




RB29H15W-110(L)
RB29H18W-110(L)
RB29H24W-110(L)




RB29H15W-110(R)
RB29H18W-110(R)
RB29H24W-110(R)



RB29H15W-102(L)
RB29H18W-102(L)
RB29H24W-102(L)




RB29H15W-102(R)
RB29H18W-102(R)
RB29H24W-102(R)




RB29H15W-F10
RB29H18W-F10
RB29H24W-F10

NOTE: When ordering units with configurations that are non-standard, drawings must accompany order.


BASE CABINETS *Sitting / Desk Height*




RB29H15W-F02
RB29H18W-F02
RB29H24W-F02




RB29H30W-040
RB29H36W-040
RB29H42W-040
RB29H48W-040




RB29H30W-050
RB29H36W-050
RB29H42W-050
RB29H48W-050




RB29H30W-044
RB29H36W-044
RB29H42W-044
RB29H48W-044




RB29H30W-006
RB29H36W-006
RB29H42W-006
RB29H48W-006



RB29H30W-060
RB29H36W-060
RB29H42W-060
RB29H48W-060




RB29H30W-200
RB29H36W-200
RB29H42W-200
RB29H48W-200




RB29H30W-210
RB29H36W-210
RB29H42W-210
RB29H48W-210


BASE CABINETS *Sitting / Desk Height*




RB29H15W-030
RB29H18W-030
RB29H24W-030




RB29H15W-B11
RB29H18W-B11
RB29H24W-B11



RB29H15W-022
RB29H18W-022
RB29H24W-022



RB29H15W-006
RB29H18W-006
RB29H24W-006



RB29H30W-100(RorL)SP
RB29H36W-100(RorL)SP
RB29H42W-100(RorL)SP
RB29H48W-100(RorL)SP



RB29H30W-220
RB29H36W-220
RB29H42W-220
RB29H48W-220




RB29H30W-110(R)SP
RB29H36W-110(R)SP
RB29H42W-110(R)SP
RB29H48W-110(R)SP




RB29H30W-110(L)SP
RB29H36W-110(L)SP
RB29H42W-110(L)SP
RB29H48W-110(L)SP

CORNER BASE CABINETS *Sitting / Desk Height*



Typical Dimensions


RC29H14W-100(L)




(SHOWN WITH OPTIONAL
LAZY SUSAN INSTALLED)

RC29H14W-100(R)


CORNER SINK BASE CABINETS *Sitting / Desk Height*




RC29H14W-10S(L)



RC29H14W-10S(R)




RC29H14W-10S(L)
Louvered



RC29H14W-10S(R)
Louvered

OPTIONAL LAZY SUSAN (For corner base cabinets)



RLS-29H
RLS-35H

WALL CASE UNITS


Wall cases are available in sizes ranging from 24", 31", 36" and 48" high and 18" wide to 48" wide. All units are 13" deep (11.75" deep I.D. w/ 10.5" D Shelves) as standard. 24" high units have one shelf. 31" high to 36" high units have two shelves. 48" high units have three shelves. Additional shelves are available. On single hinge door units, specify (R) or (L) hinge. 48" high units have 3 hinges per door.

Part Numbering System

Sample: RW30H36W-2GS 1= Wall Case 4= #Doors
 11222333-4567 2= Height 5= Material of doors*
 3= Width 6= Type of doors* *


- * Steel, Glass or no doors (open)
- ** Hinged, Sliding or Sliding/framed
- * "R" = Right hinge facing cabinet
- * "L" = Left hinge facing cabinet

7= Depth Note: 13" dp. standard,
 16" deep cabinets available -
 indicate by final suffix of 16,
 (Example: RW24H18W-00016)




OPEN SHELVING UNITS:

RW24H18W-000	RW30H18W-000	RW36H18W-000	RW48H18W-000
RW24H24W-000	RW30H24W-000	RW36H24W-000	RW48H24W-000
RW24H30W-000	RW30H30W-000	RW36H30W-000	RW48H30W-000
RW24H36W-000	RW30H36W-000	RW36H36W-000	RW48H36W-000
RW24H42W-000	RW30H42W-000	RW36H42W-000	RW48H42W-000
RW24H48W-000	RW30H48W-000	RW36H48W-000	RW48H48W-000




STEEL DOOR HINGED LEFT:

RW24H18W-1SH(L) RW30H18W-1SH(L)
 RW36H18W-1SH(L) RW48H18W-1SH(L)
 RW24H24W-1SH(L) RW30H24W-1SH(L)
 RW36H24W-1SH(L) RW48H24W-1SH(L)




STEEL DOOR HINGED RIGHT:

RW24H18W-1SH(R) RW30H18W-1SH(R)
 RW36H18W-1SH(R) RW48H18W-1SH(R)
 RW24H24W-1SH(R) RW30H24W-1SH(R)
 RW36H24W-1SH(R) RW48H24W-1SH(R)




GLASS W/ FRAME HINGED LEFT:

RW24H18W-1GH(L) RW30H18W-1GH(L)
 RW36H18W-1GH(L) RW48H18W-1GH(L)
 RW24H24W-1GH(L) RW30H24W-1GH(L)
 RW36H24W-1GH(L) RW48H24W-1GH(L)




GLASS W/ FRAME HINGED RIGHT:

RW24H18W-1GH(R) RW30H18W-1GH(R)
 RW36H18W-1GH(R) RW48H18W-1GH(R)
 RW24H24W-1GH(R) RW30H24W-1GH(R)
 RW36H24W-1GH(R) RW48H24W-1GH(R)



STEEL DOORS - HINGED:


RW24H30W-2SH RW36H30W-2SH
 RW24H36W-2SH RW36H36W-2SH
 RW24H42W-2SH RW36H42W-2SH
 RW24H48W-2SH RW36H48W-2SH
 RW30H30W-2SH RW48H30W-2SH
 RW30H36W-2SH RW48H36W-2SH
 RW30H42W-2SH RW48H42W-2SH
 RW30H48W-2SH RW48H48W-2SH



STEEL DOORS - SLIDING:


RW24H30W-2SS RW36H30W-2SS
 RW24H36W-2SS RW36H36W-2SS
 RW24H42W-2SS RW36H42W-2SS
 RW24H48W-2SS RW36H48W-2SS
 RW30H30W-2SS RW48H30W-2SS
 RW30H36W-2SS RW48H36W-2SS
 RW30H42W-2SS RW48H42W-2SS
 RW30H48W-2SS RW48H48W-2SS

WALL CASE UNITS




GLASS W/ FRAME – HINGED:

RW24H30W-2GH W36H30W-2GH
RW24H36W-2GH W36H36W-2GH
RW24H42W-2GH W36H42W-2GH
RW24H48W-2GH W36H48W-2GH
RW30H30W-2GH W48H30W-2GH
RW30H36W-2GH W48H36W-2GH
RW30H42W-2GH W48H42W-2GH
RW30H48W-2GH W48H48W-2GH




GLASS DOORS – SLIDING:

RW24H30W-2GS W36H30W-2GS
RW24H36W-2GS W36H36W-2GS
RW24H42W-2GS W36H42W-2GS
RW24H48W-2GS W36H48W-2GS
RW30H30W-2GS W48H30W-2GS
RW30H36W-2GS W48H36W-2GS
RW30H42W-2GS W48H42W-2GS
RW30H48W-2GS W48H48W-2GS




GLASS W/ FRAME – SLIDING:

RW24H30W-2GSF R W36H30W-2GSF
RW24H36W-2GSF R W36H36W-2GSF
RW24H42W-2GSF R W36H42W-2GSF
RW24H48W-2GSF R W36H48W-2GSF
RW30H30W-2GSF R W48H30W-2GSF
RW30H36W-2GSF R W48H36W-2GSF
RW30H42W-2GSF R W48H42W-2GSF
RW30H48W-2GSF R W48H48W-2GSF



PASS-THRU (OPEN ON TWO SIDES)


RWDS30H30W-00013 open-2-s
RWDS30H36W-00013 open-2-s
RWDS30H48W-00013 open-2-s



PASS-THRU (OPEN ON TWO SIDES)


RWDS30H36W-4GS13 open-2-s
RWDS30H48W-4GS13 open-2-s

CORNER WALL CASE UNITS




OPEN SHELVING:

RWC24H18W-000
RWC30H18W-000
RWC36H18W-000
RWC30H18W-00016
RWC36H18W-00016




STEEL DOOR – HINGED:

RWC24H18W-1SH(RorL)
RWC30H18W-1SH(RorL)
RWC36H18W-1SH(RorL)
RWC30H18W-1SH(RorL)16
RWC36H18W-1SH(RorL)16



GLASS DOOR W/ FRAME:

RWC24H18W-1GH(L)
RWC30H18W-1GH(L)
RWC36H18W-1GH(L)
RWC30H18W-1GH(L)16
RWC36H18W-1GH(L)16



GLASS DOOR W/ FRAME:

RWC24H18W-1GH(R)
RWC30H18W-1GH(R)
RWC36H18W-1GH(R)
RWC30H18W-1GH(R)16
RWC36H18W-1GH(R)16

NOTE: When ordering units with configurations that are non-standard, drawings must accompany order.
NOTE: 16" deep corner wall case dimension = 28-11/16"


FLOOR CASE UNITS

Floor cases are available in widths ranging from 18" to 48". All units are 84" high and are available in either 16" or 21-1/4" depths. Five shelves are standard equipment. Additional shelves are available. single hinge door units, specify (R) or (L) hinge.

Part Numbering System


Sample: RF84H18W-2SS16 1= Floor Case 4= #Doors
 11222333-456 2= Height 5= Material of
 3= Width 6= Type of doors* *
 7= Depth

- * Steel Glass or no doors (open)
- ** Hinged, Sliding or Sliding/framed
- * "R" = Right hinge facing cabinet
- * "L" = Left hinge facing cabinet




OPEN SHELVING:

RF84H18W-00016
 RF84H18W-00021
 RF84H24W-00016
 RF84H24W-00021




OPEN SHELVING:

RF84H30W-00016
 RF84H30W-00021
 RF84H36W-00016
 RF84H36W-00021
 RF84H42W-00016
 RF84H42W-00021
 RF84H48W-00016
 RF84H48W-00021




STEEL DOOR – HINGED LEFT:

RF84H18W-1SH16(L)
 RF84H18W-1SH21(L)
 RF84H24W-1SH16(L)
 RF84H24W-1SH21(L)



STEEL DOOR – HINGED RIGHT:


RF84H18W-1SH16(R)
 RF84H18W-1SH21(R)
 RF84H24W-1SH16(R)
 RF84H24W-1SH21(R)



End View (16"Dp. Unit)


NOTE: When ordering units with configurations that are non-standard, drawings must accompany order.

FLOOR CASE UNITS




GLASS W/ FRAME – HINGED LEFT:

RF84H18W-1GH16(L)
RF84H18W-1GH21(L)
RF84H24W-1GH16(L)
RF84H24W-1GH21(L)




GLASS W/ FRAME – HINGED RIGHT:

RF84H18W-1GH16(R)
RF84H18W-1GH21(R)
RF84H24W-1GH16(R)
RF84H24W-1GH21(R)




STEEL DOORS – HINGED:

RF84H30W-2SH16
RF84H30W-2SH21
RF84H36W-2SH16
RF84H36W-2SH21
RF84H42W-2SH16
RF84H42W-2SH21
RF84H48W-2SH16
RF84H48W-2SH21




STEEL DOORS – SLIDING:

RF84H30W-2SS16
RF84H30W-2SS21
RF84H36W-2SS16
RF84H36W-2SS21
RF84H42W-2SS16
RF84H42W-2SS21
RF84H48W-2SS16
RF84H48W-2SS21



GLASS W/ FRAME – SLIDING:

RF84H30W-2GSF16
RF84H30W-2GSF21
RF84H36W-2GSF16
RF84H36W-2GSF21
RF84H42W-2GSF16
RF84H42W-2GSF21
RF84H48W-2GSF16
RF84H48W-2GSF21



GLASS W/ FRAME – HINGED:

RF84H30W-2GH16
RF84H30W-2GH21
RF84H36W-2GH16
RF84H36W-2GH21
RF84H42W-2GH16
RF84H42W-2GH21
RF84H48W-2GH16
RF84H48W-2GH21

NOTE: When ordering units with configurations that are non-standard, drawings must accompany order.

FUME HOODS

QUALITY


The name RDM has long been synonymous with quality and longevity. Our newly redesigned fume hood line is a proud example of this reputation. We are backed with over 50 years of experience, combined with state of the art engineering, manufacturing and testing facilities to provide you with a top quality fume hood solution unsurpassed in quality and performance.

DESIGN

The design of our fume hood line is founded on three basic principals: safety, flexibility and aesthetics.

SAFETY

Safety should be the number one concern when choosing a fume hood. That is why our hoods are tested and certified in accordance with the Ashrae method of testing performance of laboratory fume hoods.



FLEXIBILITY

Construction of our fume hood begins with a heavy duty internal framework. This universal framework is the foundation for a variety of standard and optional features.

These features include. . .

- Chemical resistant interiors of resin-chem, PVC or stainless steel
- PVC or stainless sash guides
- Vertical, horizontal or combination sashes in a variety of materials
- Easy access for service fixture installation
- Manual or remote baffle adjustment
- Fluorescent light fixture
- Optional factory installation of service fixtures

The flexibility of our design will allow you to select standard options that will meet the requirements of your application.

AESTHETICS

Our quality, pride and craftsmanship will insure that your RDM Fume Hood will be a handsome center piece in your laboratory for years to come.



CATALOG

For a complete fume hood catalog contact your local RDM representative or contact us direct at the following address:

RDM Industrial Products, Inc.
1652 Watson Court
Milpitas, CA 95035
PH: (408) 945-8400
FX: (408) 945-8433



FUME HOODS

BYPASS & VARIABLE AIR VOLUME (VAV) FUME HOODS



PERCHLORIC ACID AND RADIO-ISOTOPE FUME HOODS


ADA FUME HOODS



WALK-IN AND SPECIAL FUME HOODS


APRONS AND DRAWER HOUSINGS

Part Numbering System: RA04H60W22D-000 1= Apron 4= Depth
 11222333444-555 2= Height 5= # Drawers
 3= Width




RA04H24W22D-000
 RA04H30W22D-000
 RA04H36W22D-000
 RA04H42W22D-000
 RA04H48W22D-000
 RA04H54W22D-000
 RA04H60W22D-000
 RA04H66W22D-000
 RA04H72W22D-000

RA04H24W28D-000
 RA04H30W28D-000
 RA04H36W28D-000
 RA04H42W28D-000
 RA04H48W28D-000
 RA04H54W28D-000
 RA04H60W28D-000
 RA04H66W28D-000
 RA04H72W28D-000




RA04H24W22D-001
 RA04H30W22D-001
 RA04H36W22D-001
 RA04H42W22D-001

RA04H24W28D-001
 RA04H30W28D-001
 RA04H36W28D-001
 RA04H42W28D-001



RA04H42W22D-002
 RA04H48W22D-002
 RA04H54W22D-002
 RA04H60W22D-002
 RA04H66W22D-002
 RA04H72W22D-002

RA04H42W28D-002
 RA04H48W28D-002
 RA04H54W28D-002
 RA04H60W28D-002
 RA04H66W28D-002
 RA04H72W28D-002




RA04H72W22D-003

RA04H72W28D-003

Note: Apron heights are 4-1/2" high.

KNEE SPACE RAILS



RKSR04H18W
 RKSR04H24W
 RKSR04H30W
 RKSR04H36W
 RKSR04H42W
 RKSR04H48W


RKSR04H52W
 RKSR04H54W
 RKSR04H60W
 RKSR04H66W
 RKSR04H72W
 RKSR04H84W
 RKSR04H96W

NOTE: When ordering units with configurations that are non-standard, drawings must accompany order.

UTILITY TABLES AND COMPONENTS

29" HIGH TABLES

35" HIGH TABLES



22" DEEP

RUT29H24W22D-000
RUT29H30W22D-000
RUT29H36W22D-000
RUT29H42W22D-000
RUT29H48W22D-000
RUT29H60W22D-000
RUT29H72W22D-000

28" DEEP


RUT29H24W28D-000
RUT29H30W28D-000
RUT29H36W28D-000
RUT29H42W28D-000
RUT29H48W28D-000
RUT29H60W28D-000
RUT29H72W28D-000

22" DEEP

RUT35H24W22D-000
RUT35H30W22D-000
RUT35H36W22D-000
RUT35H42W22D-000
RUT35H48W22D-000
RUT35H60W22D-000
RUT35H72W22D-000

28" DEEP

RUT35H24W28D-000
RUT35H30W28D-000
RUT35H36W28D-000
RUT35H42W28D-000
RUT35H48W28D-000
RUT35H60W28D-000
RUT35H72W28D-000



22" DEEP

RUT29H24W22D-001
RUT29H30W22D-001
RUT29H36W22D-001
RUT29H42W22D-001

28" DEEP


RUT29H24W28D-001
RUT29H30W28D-001
RUT29H36W28D-001
RUT29H42W28D-001

22" DEEP

RUT35H24W22D-001
RUT35H30W22D-001
RUT35H36W22D-001
RUT35H42W22D-001

28" DEEP

RUT35H24W28D-001
RUT35H30W28D-001
RUT35H36W28D-001
RUT35H42W28D-001



22" DEEP

RUT29H48W22D-002
RUT29H60W22D-002
RUT29H72W22D-002

28" DEEP


RUT29H48W28D-002
RUT29H60W28D-002
RUT29H72W28D-002

22" DEEP

RUT35H48W22D-002
RUT35H60W22D-002
RUT35H72W22D-002

28" DEEP

RUT35H48W28D-002
RUT35H60W28D-002
RUT35H72W28D-002



22" DEEP

RUT29H72W22D-003

28" DEEP

RUT29H72W28D-003

22" DEEP

RUT35H72W22D-003

28" DEEP


RUT35H72W28D-003

NOTE: UTILITY TABLES DO NOT INCLUDE COUNTER TOPS. 2"SQ. LEGS, BLACK VINYL BOOTS.

NOTE: When ordering units with configurations that are non-standard, drawings must accompany order.

ACCESSORIES

TABLE LEGS



RL29H-22D
RL35H-22D
RL29H-28D
RL35H-28D

SINGLE LEG




RSL-29H
RSL-35H

SINGLE LEG WITH PLATE




RSL-29HP
RSL-35HP

END PEDESTAL




REP29H02W-22D
REP35H02W-22D
REP29H02W-28D
REP35H02W-28D

WALL CASE (UPPER CAB) FRONT FILLER



Specify Width (W)
RWCFF24H-W
RWCFF30H-W
RWCFF36H-W
RWCFF48H-W

WALL CASE CORNER FILLER




Specify Depth (D)
13" or 16"


RWCCF24H-D
RWCCF30H-D
RWCCF36H-D
RWCCF48H-D

Angle is 1"

WALL - FLOOR CASE SHELF



REAGENT SHELF UPRIGHT



RRSU18H12D2W
(18"Hx12"Dx2"W)

Note: Special heights & options available upon request. Drawings and dimensions must accompany order.

NOTE: Use with RSF Fascia Panels for Lateral Top Support (Not Shown) - Mounts between Uprights.

SLOPE TOP WITH END CAPS


PART NUMBERS

WALL

RWS18W13D
RWS24W13D
RWS30W13D
RWS36W13D
RWS42W13D
RWS48W13D
RWS18W16D
RWS24W16D
RWS30W16D
RWS36W16D
RWS42W16D
RWS48W16D

FLOOR

RFS18W16D
RFS24W16D
RFS30W16D
RFS36W16D
RFS42W16D
RFS48W16D
RFS18W21D
RFS24W21D
RFS30W21D
RFS36W21D
RFS42W21D
RFS48W21D




RST13D-18W RST21D-18W
RST13D-24W RST21D-24W
RST13D-30W RST21D-30W
RST13D-36W RST21D-36W
RST13D-42W RST21D-42W
RST13D-48W RST21D-48W
RST16D-18W
RST16D-24W
RST16D-30W
RST16D-36W
RST16D-42W
RST16D-48W

NOTE: When ordering units with configurations that are non-standard, drawings must accompany order.

ACCESSORIES

KNEE SPACE PANELS



For use in sit down situation with apron

24" HIGH

RKSP24H24W
RKSP24H30W
RKSP24H36W
RKSP24H42W
RKSP24H48W
RKSP24H54W
RKSP24H60W
RKSP24H66W
RKSP24H72W

For use in standing situation with apron, or sit down without apron

29" HIGH

RKSP29H24W
RKSP29H30W
RKSP29H36W
RKSP29H42W
RKSP29H48W
RKSP29H54W
RKSP29H60W
RKSP29H66W
RKSP29H72W

For use in standing situation without apron

35" HIGH

RKSP35H24W
RKSP35H30W
RKSP35H36W
RKSP35H42W
RKSP35H48W
RKSP35H54W
RKSP35H60W
RKSP35H66W
RKSP35H72W

FILL PANELS

All base cabinet fillers are 35" high. Fillers for 29" high (low base) cabinets will need to be trimmed by the installer. Other sizes available upon request.




C-FILLER

RCF35-8W
RCF35-14W



L-FILLER

RLF35-1W
RLF35-7W
RLF35-11W



C-FRONT FILLER

RFFC35-8W
RFFC35-14W



CORNER FILLER

RICF35H-1W
RICF35H-2W
RICF35H-4W



L-FRONT FILLER

RFFL35H-8W
RFFL35H-14W



FLOOR CASE FILLER

RFFL84H-2W
RFFL84H-8W



TOE SPACE FILLER

RTSF-00 (WIDTH)

* Last number(s) in filler designation indicates width of filler assembly.



Represented By:



RDM LABORATORY FURNITURE SOLUTIONS

2020 Traders Ridge Drive

Conroe, Texas 77301

Telephone: (877) 777-9130

FAX: (936) 777-8707

Web page: www.rdm-ind.com

E-mail: lab-sales@rdm-ind.com

© RDM Industrial Products, Inc.

2025 Printed in the U.S.A.

QUALITY • SERVICE • TIMING